---

Gateway)HV<sup>®</sup>  
Hemostasis Valve Adapter

---

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	5
Anwendungsgebiete	6
Ενδείξεις χρήσης	7
Felhasználási javallatok	8
Indicazioni per l'uso	9
Bruksindikasjoner	10
Indicações de utilização	11
Indicaciones de uso	12
Indikationer för användning	13
Kullanım Endikasyonları	14

Gateway<sup>®</sup> HV<sup>™</sup>  
Hemostasis Valve Adapter

### Hemostasis Valve Adapter

**Indications for Use:**

The Hemostasis Valve Adapter is indicated to minimize blood loss during the introduction of catheters and other intravascular devices into the vasculature.

**Intended Patient Population and Group:**

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

**Device Description:**

The Hemostasis Valve Adapter consists of a hub with a locking collar, luer fitting, hemostasis valve, and side arm tubing with stopcock.

**Contraindications:**

None known

**Potential Complications:**

The potential complications related to the use of the valve adapter include, but are not limited to the following: Blood Loss, Air embolism.

**Precautions:**

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

**USA CAUTION**

**Federal (USA) law restricts this device to sale by or on the order of a physician.**

**Cautions:**

- Read instructions prior to use.
- Use of this device should only be performed by physicians thoroughly trained in this procedure.
- Do not expose to organic solvents, eg. Alcohol. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Damage to the valve may occur if catheter or other device is not centered in the valve when inserted, left in the valve for extended periods or withdrawn rapidly from the device.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

**Warnings:**

- Do not alter this device in any way.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not resterilize.
- The Hemostasis Valve Adapter may be used with catheters or other devices that are 8F or smaller.
- Do not exceed 30 psi internal pressure.

**USE STERILE TECHNIQUE.** A suggested procedure:

1. Peel open package and place contents on sterile field. Inspect device and accessories for defects. Do not use any defective devices.
2. Follow standard hospital practice for vascular access.
3. The Hemostasis Valve Adapter may be connected to the catheter or device prior to insertion into the vascular or after the catheter is in place.
4. Follow manufacturer's recommendation is for the catheter or device being introduced through the Hemostasis Valve Adapter.

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