Needle

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	6
Anwendungsgebiete	8
Ενδείξεις χρήσης	10
Felhasználási javallatok	12
Indicazioni per l'uso	13
Bruksindikasjoner	15
Indicações de utilização	16
Indicaciones de uso	17
Indikationer för användning	18
Kullanım Endikasyonları	19

Needle

Indications for Use:

These needles are used for the percutaneous introduction of guidewires.

Intended Patient Population and Group:

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

Device Description:

Introducer Needles are composed of two components: a stainless steel cannula with an over molded hub. The introducer needle provides an access path into the vasculature.

Contraindications:

Use of the needle is contraindicated if the patient has a known or suspected obstruction in the vessel.

Potential Complications:

The potential complications related to the use of the needle include, but are not limited to the following: Air embolism, device dislodgement, needle cannula breakage, hematoma, hemothorax, pneumothorax, thrombosis, trauma to vessels, sepsis.

Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

Cautions:

- This procedure should only be performed by professionals thoroughly trained in this procedure.
- If resistance is met when advancing or withdrawing the guidewire through the needle, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

Warnings:

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

USE STERILE TECHNIQUE. A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated puncture site as desired.
3. Insert appropriate size Needle into vessel – observe for flash back.
4. Follow standard hospital practice for placement of guidewire.

Needle Disposal: The needle should be disposed of in accordance with policies and procedures for contaminated sharps disposal. Single use only.

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