Sharps Safety Guidewire Introducer Needle
SP-Pro™ GWI

Indications for Use:

These needles are used for the percutaneous introduction of guidewires.

Intended Patient Population and Group:

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

Device Description:

The SP-Pro™ GWI Guide Wire Introducer is a single use device for introduction of guidewires into the vascular system. The device is designed with a tip shield and sheath mechanism integrated with the introducer hub to aid in the prevention of sharps injuries. The shield covers the sharp tip of the device by actively disconnecting it from the handle and pulling the introducer out of the body, off the guidewire and back into the shield. After activation the shield covers the introducer tip.

Contraindications:

Use of the needle is contraindicated if the patient has a known or suspected obstruction in the vessel.

Potential Complications:

The potential complications related to the use of the needle include, but are not limited to the following: Air embolism, device dislodgement, needle cannula breakage, hematoma, hemothorax, pneumothorax, thrombosis, trauma to vessels, sepsis.

Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

Cautions:

- This procedure should only be performed by professionals thoroughly trained in this procedure.
- If resistance is met when advancing or withdrawing the guidewire through the needle, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

Warnings:

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

USE STERILE TECHNIQUE. A suggested procedure:

Preparation for Use

1. Open package and place needle on sterile field.
2. Grip the clear safety handle and remove needle guard.
3. Follow standard hospital practice for veni-puncture and placement of guide-wire.
4. Remove SP-Pro™ needle from vessel and guide-wire. Keep hands behind needle at all times prior to deployment of safety feature.

Deploying the Sharps Injury Prevention Feature

Step 5

Step 6

Step 7

5. While holding the safety handle (luer toward user), rotate the needle flashback chamber counter-clockwise to unlock the safety feature and pull back.
6. Withdraw the needle until the needle tip disappears into the safety handle and locks securely into needle handle as indicated by an audible "click" and feel. Hearing loss and dexterity may hamper these signs of activation. Follow the next step carefully!
7. Deployment of Sharps Injury Prevention Feature is confirmed by gently pushing the needle flashback chamber forward to verify that the needle is securely locked in the needle handle. If needle is not locked into handle repeat locking procedure. If for any reason safety mechanism does not deploy use extreme caution and proceed to step 8. Please report any occurrence of non-deployment to Galt Medical Corp.

Needle Disposal

8. The needle should be disposed of in accordance with policies and procedures for contaminated sharps disposal. Do not attempt to override or defeat the locking safety mechanism. Single Use Only.

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**Sharps Safety Guidewire Introducer Needle
SP-Pro™ GWI**

Outer Diameter
 udvendig diameter
 Buitendiameter
 Diamètre extérieur
 Außendurchmesser
 Εξωτερική διάμετρος
 Külső átmérő
 Diametro esterno
 ytre Diameter
 Diámetro exterior
 Diámetro exterior
 ytterdiameter
 Dis Cap

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Emergo Europe
 Westervoortsedijk 60
 6827 AT Arnhem
 The Netherlands

CH REP

MedEnvoy Switzerland
 Gotthardstrasse 28
 6302 Zug
 Switzerland

SP-Pro™ GWI is manufactured by NeoMedical Inc., and protected by the following patents:
 U.S. Patents 8,066,678 & 7,717,888; Europe EP 1430921; Japan JP 4355567

Rx ONLY

Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner)

STERILE EO

Ethylene Oxide Sterilized

Keep Dry

REF

Part Number

Instructions for Use

LOT

Lot Number

Single Use Only

Caution

Manufacturing Date

Quantity

Not Made With Natural Rubber Latex

Non-Pyrogenic

Use By Date

MD

Medical Device

Do Not Re-Sterilize

Store In A Dark Place

Importer

MR Unsafe

Manufacturer

EU Authorized Representative

Swiss Authorized Representative

NiTi

Nitinol

NiTi [Au]

Nitinol Gold

NiTi [Pt]

Nitinol Platinum

Single Sterile Barrier system

Single Sterile Barrier System with Protective Packaging Outside