---

**Trocar Needle**

---

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	6
Anwendungsgebiete	8
Ενδείξεις χρήσης	10
Felhasználási javallatok	12
Indicazioni per l'uso	13
Bruksindikasjoner	15
Indicações de utilização	16
Indicaciones de uso	18
Indikationer för användning	20
Kullanım Endikasyonları	21

---

## Trocar Needle

---

### Indications for Use:

These needles are used for the percutaneous introduction of guidewires.

### Intended Patient Population and Group:

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

### Device Description:

Introducer Needles are composed of two components: a stainless steel cannula with an over molded hub. The introducer needle provides an access path into the vasculature.

### Contraindications:

Use of the needle is contraindicated if the patient has a known or suspected obstruction in the vessel.

### Potential Complications:

The potential complications related to the use of the needle include, but are not limited to the following: Air embolism, device dislodgement, needle cannula breakage, hematoma, hemothorax, pneumothorax, thrombosis, trauma to vessels, sepsis.

### Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

### USA CAUTION

**Federal (USA) law restricts this device to sale by or on the order of a physician.**

### Cautions:

- This procedure should only be performed by professionals thoroughly trained in this procedure.
- If resistance is met when advancing or withdrawing the guidewire through the needle, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State

### Warnings:

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

### USE STERILE TECHNIQUE. A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated puncture site per standard practice.
3. Insert appropriate size Needle into vessel – observe for flash back.
4. Follow standard clinical practice for placement of trocar.

**Needle Disposal:** The needle should be disposed of in accordance with policies and procedures for contaminated sharps disposal. Single use only.

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