---

**Elite HV® Radial Artery Access Introducer Kit**

---

Indications for Use	2
Indikationer for anvendelse	4
Gebruiksindicaties	6
Indications d'emploi	8
Indikationen	10
Ενδείξεις χρήσης	12
Felhasználási javallatok	14
Indicazioni per l'uso	16
Bruksindikasjoner	18
Indicações de utilização	20
Indicaciones de uso	22
Indikationer för användning	24
Kullanım Endikasyonları	26

---

## Elite HV® Radial Artery Access Introducer Kit

---

### Indications for Use:

The Galt Medical Corp Elite HV® Radial is indicated to facilitate placing a catheter through the skin into a vein or artery including but not limited to the radial artery

### Intended Patient Population and Group:

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

### Device Description:

The Radial Artery Access Introducer device assembly consists of an outer introducer sheath with an integrated hemostasis valve, with an optional side port, and a dilator.

Kit configurations may include the following devices:

**Guidewire:** There are two configurations; Mandrel or Fully Coiled. Mandrel guidewire construction uses a tapered core wire which is soldered or welded to a coil at the distal end. Fully Coiled guidewire construction uses a tightly wound coil which surrounds a tapered core and a safety ribbon.

**Needle:** Introducer needles are composed of two components: a stainless steel cannula with an over molded hub. The introducer needle provides an access path into the vasculature.

### Contraindications:

Use of the introducer is contraindicated if the patient has a known or suspected obstruction in the vessel, including abnormal Allen's test, radial pulse, insufficient dual arterial supply, or inadequate circulation to the hand. There is increased risk of pneumothorax for the patient who has severe chronic lung disease.

### Potential Complications:

The potential complications related to the use of the introducer include, but are not limited to the following: Air embolism, arterial spasm, arterial occlusion, arterial injury/perforation/dissection, device dislodgement, hematoma, hemorrhage, hemothorax, sepsis pneumothorax, thrombosis, peripheral nerve damage.

### Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

### USA CAUTION

**Federal (USA) law restricts this device to sale by or on the order of a physician.**

### Cautions:

- This procedure should only be performed by physicians thoroughly trained in the technique of percutaneous vessel entry and familiar with the complications of angiographic or percutaneous procedures.
- Before use, make sure the sheath (Fr.) and dilator size are appropriate for the access vessel and the system to be used.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- If resistance is met when advancing or withdrawing the guidewire or the introducer, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Use fluoroscopic guidance for placement longer sheaths and guidewires, used to navigate tortuous vessels or to position a catheter tip at a specific location.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Do not expose to organic solvents, eg. alcohol. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not attempt to use a guidewire larger than maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Apply appropriate anticoagulant therapy to the patient.
- Do not leave the introducer in place for extended periods of time without a catheter or an obturator to support the cannula wall.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State
- Damage to the valve assembly may occur under the following circumstances:
  - Obturator or catheter in valve system for extended periods.
  - Inner catheter is withdrawn too rapidly.

### Warnings:

- Do not alter this device in any way.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.
- Prior to beginning Radial Artery access, an assessment such as the Allen's test should be performed to access the presence/adequacy of dual arterial circulation to the hand.

### USE STERILE TECHNIQUE. A suggested procedure:

1. Peel open package and place contents on sterile field. Inspect Radial Artery Access Introducer and accessories for defects. Do not use any defective devices.
2. To remove air, flush the dilator, Radial Artery Access Introducer and sideport with normal saline solution.
3. Prep skin and drape in area of anticipated puncture site as desired.
4. Insert introducer needle into vessel. The needle position should be verified by observing blood return.
5. The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person.
6. Aspirate the puncture needle using a syringe.
7. Remove the syringe and insert soft tip of the guidewire through the introducer needle into the vessel. Advance guidewire to required depth. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire location is suggested.
8. Hold guidewire in place and remove introducer needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.

**CAUTION:** Do not allow Guidewire to advance totally into patient.

9. Enlarge cutaneous puncture site with cutting edge of scalpel, cutting away from guidewire.

**PRECAUTION:** Do not cut guidewire.

10. Assemble the introducer set by carefully inserting the dilator completely into sheath introducer. Firmly push the snap fit ring on the dilator into the sheath valve cap.
11. When using an introducer with a sideport, follow standard hospital practice for using a continuous drip of normal saline solution through the sideport while the hemostasis introducer is in the vessel.
12. Introducers with the hydrophilic coating should be hydrated by immersing introducer in normal saline solution.
13. While holding the radial artery access introducer set close to the skin, advance the dilator and sheath together with a twisting motion over the guidewire and into the vessel. Fluoroscopic observation may be

LAB-333-00, Rev. G (02/20)

advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.

14. To detach the dilator from the sheath cap, push the dilator hub to one side until it becomes detached. Remove the vessel dilator and guidewire, leaving the sheath as a conduit into the vessel.
15. Introduce the selected catheter or other device into the sheath using the instructions provided by the manufacturer of the catheter or other device, and standard hospital practice.
16. To change catheters, slowly withdraw the catheter from the vessel using instructions provided by the manufacturer of the catheter or other device, and standard hospital practice and repeat the insertion procedure.

**CAUTION:** When removing the catheter, aspirate via the sideport extension to collect fibrin that may have been deposited at the tip of the sheath.

**OBTURATOR, insertion and withdrawal:**

1. For occlusion of sheath, use obturator of same size as sheath.
2. For flushing and infusion, use an obturator one french size smaller than the designated sheath size.
3. While holding the sheath in place, slide the obturator through the sheath valve and firmly push the snap fit ring into the sheath valve cap.
4. Attach a flushing line to the sheath sideport extension and flush.
5. Establish an infusion drip via the sideport extension per hospital protocol.
6. While holding the sheath in place, place thumb and index finger between the sheath and obturator hubs, squeeze firmly and withdraw.

**CAUTION:** When removing the obturator, aspirate via the sideport extension to collect any fibrin that may have deposited at the tip of the sheath.

**SHEATH REMOVAL:**

1. The sheath may be removed when clinically indicated. Compression on the vessel, above the puncture site, should be started as the sheath is slowly removed.

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