

**Theragenics**  
**POLYPROPYLENE**  
**NONABSORBABLE SURGICAL**  
**SUTURES, U.S.P.**

**DESCRIPTION**

Polypropylene is a monofilament, nonabsorbable, sterile surgical suture composed of a strand of polypropylene, a synthetic linear Polyolefin. Polypropylene sutures are undyed or dyed blue (Pthalocyaninato [2-] Copper).

Polypropylene sutures meet all the requirements established by the United States Pharmacopeia (U.S.P.) for nonabsorbable surgical sutures.

**INDICATIONS**

Polypropylene is indicated for use in general soft tissue approximation/ligation.

**ACTIONS**

Polypropylene sutures elicit a minimal inflammation reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Polypropylene suture is not absorbed, nor is any significant change in tensile strength retention known to occur in vivo.

**CONTRAINDICATIONS**

The use of this suture is contraindicated on patients with known sensitivities or allergies to its components.

**WARNINGS**

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing polypropylene suture for wound closure, as a risk of wound dehiscence may vary with the site of application and the suture material used.

Do not use if package is open or damaged or if the expiration date has been exceeded.

Discard open unused suture. Do not resterilize; reesterilization may alter the physical properties of this suture, which will result in an adverse patient reaction.

Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in a "sharps" container.

Store in a cool dry environment.

As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation.

Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds.

**PRECAUTIONS**

In handling this or any other surgical suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Adequate knot security requires the accepted surgical technique of flat, squared ties with additional throws as warranted by surgical circumstance and experience of the surgeon. The use of additional throws may be particularly appropriate when tying monofilaments.

**ADVERSE REACTIONS**

Adverse effects associated with the use of this device include: wound dehiscence, calculus formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, infected wounds, minimal acute inflammatory tissue reaction, and transitory local irritation.

Discuss the potential for allergic reaction in patients that are known to be sensitive to polypropylene suture.


**HOW SUPPLIED**

Polypropylene sutures are available in various USP sizes. Polypropylene sutures are supplied in a wide range of lengths affixed to a diverse assortment of needle types.

**CAUTION**

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

**SYMBOL DEFINITIONS**

-  Lot Number
-  Expiration Date
-  Do Not Reuse
-  Do Not Resterilize
-  See Instructions For Use
-  Sterilized By Ethylene Oxide
-  Keep away from sunlight and heat
-  Do not use if package is damaged
-  Keep dry
-  Manufacturer

 **CP Medical Inc.**  
 1775 Corporate Drive, Suite 150  
 Norcross, GA 30093 USA

**Theragenics Corporation®**  
 5203 Bristol Industrial Way  
 Buford, GA 30518 USA  
 877-444-7333  
 www.theragenics.com

 PUB511 R01  
 Effective Date 08/23/2022