---

**DILATOR**

---

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	6
Anwendungsgebiete	8
Ενδείξεις χρήσης	10
Felhasználási javallatok	12
Indicazioni per l'uso	13
Bruksindikasjoner	14
Indicações de utilização	15
Indicaciones de uso	16
Indikationer för användning	17
Kullanım Endikasyonları	18

---

## DILATOR

---

### Indications for Use:

These dilators are used for the percutaneous introduction of guidewires into the peripheral vasculature.

### Intended Patient Population and Group:

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

### Device Description:

The dilator introducer consists of a tipped radiopaque tube with a molded hub that accepts up to a .038 inch guidewires

### Contraindications:

Use of the dilator is contraindicated if the patient has a known or suspected obstruction in the vessel. There is increased risk of pneumothorax for the patient who has severe chronic lung disease. Poor healing may result in the patient who has had irradiation to the anterior chest.

### Potential Complications:

The potential complications related to the use of the dilator include, but are not limited to the following: Air embolism, device embolism, device dislodgement, pneumothorax, vein thrombosis, hematoma formation, hemothorax, vessel erosion, trauma to vessels, sepsis.

### Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

### USA CAUTION

**Federal (USA) law restricts this device to sale by or on the order of a physician.**

### Cautions:

- This procedure should only be performed by professionals thoroughly trained in this procedure.
- If resistance is met when advancing or withdrawing the guidewire through the needle, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

### Warnings:

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

### USE STERILE TECHNIQUE. A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated puncture site as desired.
3. Distend the vessel of choice following standard hospital practice for venipuncture. The vein will be much easier to locate if the patient is well hydrated.
4. Insert needle into vessel. The needle position should be verified by observing venous blood return.
5. The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person.
6. Aspirate the puncture needle using the syringe.
7. Remove the syringe and insert soft tip of the guidewire through the introducer needle into the vessel. Advance guidewire to required depth. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire position is suggested.
8. Hold guidewire in place and remove introducer needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.  
**Caution: If the guidewire must be withdrawn while the needle is inserted, remove both the needle and wire as a unit to prevent the needle from damaging or shearing the guidewire.**
9. Thread the dilator over the guidewire.
10. Advance the dilator with a twisting motion over the guidewire and into the vessel, to dilate tissue in vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.
11. Slowly withdraw and remove the dilator, while holding the guidewire in position and insert catheter or over guidewire.
12. Remove guidewire leaving catheter in position.

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