---

**GalfVTI® Valved Tearaway Sheath Introducer**

---

Indications for Use	2
Indikationer for anvendelse	4
Gebruiksindicaties	6
Indications d'emploi	8
Anwendungsgebiete	10
Ενδείξεις χρήσης	12
Felhasználási javallatok	14
Indicazioni per l'uso	16
Bruksindikasjoner	18
Indicações de utilização	20
Indicaciones de uso	22
Indikationer för användning	24
Kullanım Endikasyonları	26

---

## GaltVTI® Valved Tearaway Sheath Introducer

---

### Indications for Use:

The introducer system is intended for use in percutaneous procedures to introduce catheters and other intravascular devices into the peripheral vasculature.

### Intended Patient Population and Group:

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

### Device Description:

The GaltVTI® Valved Tearaway Introducer device consists of an outer tearaway sheath with an integrated hemostasis valve, with an optional side port, and a dilator. The tearaway sheath has a winged hub to initiate the tear.

### Contraindications:

Use of the introducer is contraindicated if the patient has a known or suspected obstruction in the vessel. There is increased risk of pneumothorax for the patient who has severe chronic lung disease. Poor healing may result in the patient who has had irradiation to the anterior chest.

### Potential Complications:

The potential complications related to the use of the introducer include, but are not limited to the following: Air embolism, device embolism, wound infection, sepsis, pyrogenic response, intimal tear, artery puncture, Hematoma, hemothorax, Vessel Trauma, Vessel erosion, pneumothorax, thrombosis, device dislodgment, blood loss, patient injury/death.

### Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

### USA CAUTION

**Federal (USA) law restricts this device to sale by or on the order of a physician.**

### Cautions:

- This procedure should only be performed by physicians thoroughly trained in this procedure.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- If resistance is met when advancing or withdrawing the guidewire or the introducer, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Do not use forceps to break the handle and/or to peel the sheath as this may damage the sheath and cause premature withdrawal of the sheath from the patient.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Symmetrical peeling of the sheath is critical.
- Dilators and catheters should not remain in the sheath for prolonged periods of time to maintain the performance of the valve.
- Dilators and other device should be inserted into the center of the valve to prevent damage to the valve.
- Dilators and catheters should be removed slowly from the sheath to prevent damage to the valve.
- For sideport model, remove the dilator prior to using the sideport. Use sideport only when injecting or aspirating through sheath.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

### Warnings:

- Do not alter this device in any way.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.
- Do not use sideport prior to removing all air.
- Do not use with power injection systems.

### USE STERILE TECHNIQUE. A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated veni-puncture as desired.
3. Distend the vessel of choice following standard hospital practice for venipuncture. The vein will be much easier to locate if the patient is well hydrated.
4. Insert needle into vessel. The needle position should be verified by observing venous blood return.
5. The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person.
6. Aspirate the puncture needle using the syringe.
7. Remove the syringe and insert soft tip of the guidewire through the introducer needle into the vessel. Advance guidewire to required depth. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire position is suggested.
8. Hold guidewire in place and remove introducer needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.
9. Remove the dilator from the shipping position and insert through the proximal end of the sheath. Ensure that the dilator is inserted into the center of the valve to prevent damage to the valve. Insert all the way in rotating clockwise when almost fully inserted until you feel resistance as the dilator hub locks against the sheath hub.
10. Prime the introducer prior to assembly over the guidewire. For sideport model, flush through the sideport. For non sideport model, flush through dilator/sheath assembly. NOTE: never stick a needle directly into the sheath valve without a dilator in place.
11. Thread the dilator/sheath assembly over the guidewire.
12. Advance the dilator and sheath together with a twisting motion over the guidewire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.
13. Remove the vessel dilator by rotating counterclockwise a quarter turn while removing the guidewire, leaving the sheath as a conduit into the vessel. The GaltVTI® will minimize blood loss and air aspiration through the valve.
14. Advance the catheter or device through the center of the valve and sheath and into the vessel.

**Note:** If desired, lubricate the device with a light coating of appropriate sterile lubricant to facilitate passage through the introducer sheath. Apply moderate pressure until the device tip pushes through the distal end of the sheath. Some resistance may be felt when advancing a device having tines or a flange tip.

15. On sideport model, flush with saline prior to peeling sheath to minimize backbleeding. When the lead, catheter or other device is properly located, peel the sheath away from the device by using the ends of the sheath as handles and pulling them so that the sheath separates longitudinally and breaks at the distal end.

**Caution:** Do not attempt to withdraw a tined or flanged device through the introducer sheath. If the device must be withdrawn, remove the sheath first and then gently pull the device out of the vessel. Repeat the insertion procedure with a new introducer.

**Caution:** When the guidewire is in a vessel, do not advance the movable core if the tip is in a curved shape. Never twist or force the core because excessive force may cause it to penetrate the coil and damage the vessel.

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