Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	5
Anwendungsgebiete	6
Ενδείξεις χρήσης	7
Felhasználási javallatok	8
Indicazioni per l'uso	9
Bruksindikasjoner	10
Indicações de utilização	11
Indicaciones de uso	12
Indikationer för användning	13
Kullanım Endikasyonları	14

---

## GUIDEWIRES

---

**Indications for Use:**

Guidewires are intended for use in percutaneous procedures, to introduce and position catheters and other interventional devices within the coronary and peripheral vasculature.

**Intended Patient Population and Group:**

For use on both adult and adolescent groups. The specific target population is dependent on the delivered device.

**Device Description:**

Guidewires come in two configurations; Mandrel or Fully Coiled. Mandrel guidewire construction uses a tapered core wire which is soldered or welded to a coil at the distal end. Fully Coiled guidewire construction uses a tightly wound coil which surrounds a tapered core and a safety ribbon.

**Contraindications:**

Use of the guidewire is contraindicated if the patient has a known or suspected obstruction in the vasculature.

**Potential Complications:**

The potential complications related to the use of the guidewire include, but are not limited to the following: Air embolism, wound infection, hematoma at the puncture site, perforation of the vessel wall vessel trauma, thrombosis, device embolism, cardiac tamponade, cardiac arrhythmia, and guidewire break or unravel.

**Precautions:**

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

**USA CAUTION**

**Federal (USA) law restricts this device to sale by or on the order of a physician.**

**Cautions:**

- This procedure should only be performed by physicians thoroughly trained in this procedure.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- If resistance is met when advancing or withdrawing the guidewire, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

**Warnings:**

- Do not alter this device in any way.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear, unravel, or PTFE coating may scrape off.
- Do not resterilize.

**USE STERILE TECHNIQUE, A suggested procedure:**

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated veni-puncture as desired.
3. Locate vessel using small gauge needle and syringe (not supplied).
4. Insert appropriate gauge Introducer Needle into vessel - observe for flash back.
5. Flush the guidewire with sterile heparinized normal saline or a similar isotonic solution. Insert desired flexible tip of the guidewire through the introducer needle into vessel. Slowly advance guidewire to required depth. Advancement of "J" tip may require a gentle rotating motion.  
**Caution:** At no time should guidewire be advanced or withdrawn when resistance is met. Determine cause of resistance before proceeding.
6. Hold guidewire in place and remove Introducer Needle.  
**Caution:** Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.
7. With the guidewire in place follow the instructions for use provided by the manufacturer of the introducer, catheter or other device, as well as the medical facility's standard procedure for these types of devices.
8. To lengthen the flexible segment of a movable core guidewire, gently withdraw the core. To shorten the flexible segment, push the core in carefully.  
**Caution:** When the guidewire is in a vessel, do not advance the movable core if the tip is in a curved shape. Never twist or force the core because excessive force may cause it to penetrate the coil and damage the vessel.

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
**GUIDEWIRES**

<b>Rx ONLY</b>				
Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner)	Do Not Use if Package Is Damaged	Single Use Only	Not Made With Natural Rubber Latex	Do Not Re-Sterilize
				
Ethylene Oxide Sterilized	Instructions for Use	Caution	Non-Pyrogenic	Store In A Dark Place
				
Keep Dry	Part Number	Lot Number	Use By Date	Importer
				
Manufacturer	EU Authorized Representative	Swiss Authorized Representative	Medical Device	MR Unsafe
<b>NiTi</b>	<b>NiTi [Au]</b>	<b>NiTi [Pt]</b>		
Nitinol	Nitinol Gold	Nitinol Platinum	Single Sterile Barrier system	Single Sterile Barrier System with Protective Packaging Outside