MicroSlide® Tearaway Sheath Introducer

Indications for Use	2
Indikationer for anvendelse	4
Gebruiksindicaties	6
Indications d'emploi	8
Anwendungsgebiete	10
Ενδείξεις χρήσης	12
Felhasználási javallatok	14
Indicazioni per l'uso	16
Bruksindikasjoner	18
Indicações de utilização	20
Indicaciones de uso	22
Indikationer för användning	24
Kullanım Endikasyonları	26

MicroSlide® Tearaway Sheath Introducer

Indications for Use:

These introducers are used for the procedures to introduce catheters and other intravascular devices into the coronary and peripheral vasculature of adult and pediatric patients of all ages.

Device Description:

The MicroSlide® device consists of an outer tearaway sheath and a dilator. The tear-away sheath has a winged hub to initiate the tear.

Kit configuration includes:

Guidewire: There are two configurations; Mandrel or Fully Coiled. Mandrel guidewire construction uses a tapered core wire which is soldered or welded to a coil at the distal end. Fully Coiled guidewire construction uses a tightly wound coil which surrounds a tapered core and a safety ribbon.

Needle: Introducer needles are composed of two components: a stainless steel cannula with an over molded hub. The introducer needle provides an access path into the vasculature.

Contraindications:

Use of the introducer is contraindicated if the patient has a known or suspected obstruction in the vessel. There is increased risk of pneumothorax for the patient who has severe chronic lung disease. Poor healing may result in the patient who has had irradiation to the anterior chest.

Potential Complications:

The potential complications related to the use of the introducer include, but are not limited to the following: Air embolism, device embolism, device dislodgement, pneumothorax, vein thrombosis hematoma formation, hemothorax, vessel erosion trauma to vessels, sepsis.

Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

Cautions:

- This procedure should only be performed by clinicians thoroughly trained in this procedure.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- If resistance is met when advancing or withdrawing the guidewire or the introducer, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Do not use forceps to break the handle and/or to peel the sheath as this may damage the sheath and cause premature withdrawal of the sheath from the patient.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and clinician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Symmetrical peeling of the sheath is critical.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

Warnings:

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

Pediatric Precautions:

When Used in the pediatric population, adhere to the following precautions:

Additional Contraindications to consider:

- Presence of Bacteremia or septicemia (or is suspected), Vessel size is insufficient to accommodate the size of the introducer, Anatomical irregularities could interfere with introducer placement, Patient is known or suspected to be allergic to materials contained in the device, Disruption of local tissue at access site, (revision of current) Patients who have had or anticipated having radiation to site.

Additional Complications to consider:

- Bleeding, Nerve injury, Hematoma, Vessel laceration/perforation, Thromboembolism.

Additional Cautions to consider:

- Do not advance the guidewire past the axillary region.

USE STERILE TECHNIQUE. A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated venipuncture as desired.
3. Distend the vessel of choice following standard hospital practice for venipuncture. The vein will be much easier to locate if the patient is well hydrated.
4. Insert needle into vessel. The needle position should be verified by observing venous blood return.
5. The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person.
6. Aspirate the puncture needle using the syringe.
7. Remove the syringe and insert soft tip of the guidewire through the introducer needle into the vessel. Advance guidewire to required depth. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire tip is suggested.
8. Hold guidewire in place and remove introducer needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.
9. Thread the dilator/sheath assembly over the guidewire.
10. Advance the dilator and sheath together with a twisting motion over the guidewire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.
11. Remove the vessel dilator and guidewire, leaving the sheath as a conduit into the vessel. Immediately place a finger over the remaining sheath orifice to prevent excessive bleeding or possible air aspiration.
12. Advance the catheter or device through the sheath and into the vessel.

Note: If desired, lubricate the device with a light coating of appropriate sterile lubricant to facilitate passage through the introducer sheath. Apply moderate pressure until the device tip pushes through the distal end of the sheath. Some resistance may be felt when advancing a device having tines or a flange tip.

13. When the device, catheter or other device is properly located, peel the sheath away from the device by using the ends of the sheath as handles and pulling them so that the sheath separates longitudinally and breaks at the distal end.

Caution: Do not attempt to withdraw a tined or flanged device through the introducer sheath. If the device must be withdrawn, remove the sheath first and then gently pull the lead out of the vessel. Repeat the insertion procedure with a new introducer.

Caution: When the guidewire is in a vessel, do not advance the movable core if the tip is in a curved shape. Never twist or force the core because excessive force may cause it to penetrate the coil and damage the vessel.

MicroSlide® Tearaway-sheathintroducer

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
Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

<p>Rx ONLY</p> <p>Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner)</p>	 <p>Do Not Use if Package Is Damaged</p>	 <p>Single Use Only</p>	 <p>Not Made With Natural Rubber Latex</p>	 <p>Do Not Re-Sterilize</p>
 <p>Ethylene Oxide Sterilized</p>	 <p>Instructions for Use</p>	 <p>Caution</p>	 <p>Non-Pyrogenic</p>	 <p>Store In A Dark Place</p>
 <p>Keep Dry</p>	 <p>Part Number</p>	 <p>Lot Number</p>	 <p>Manufacturing Date</p>	 <p>Use By Date</p>
 <p>Manufacturer</p>	 <p>EU Authorized Representative</p>	 <p>Swiss Authorized Representative</p>	 <p>Quantity</p>	 <p>Medical Device</p>
<p>NiTi Nitinol</p>	<p>NiTi [Au] Nitinol Gold</p>	<p>NiTi [Pt] Nitinol Platinum</p>	 <p>Single Sterile Barrier system</p>	 <p>Single Sterile Barrier System with Protective Packaging Outside</p>

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Galt Medical, Corp.