



White paper

Focaltherapy for localised prostate cancer: systematic review of efficacy and functional outcomes from partial gland ablations with Low Dose Rate Brachytherapy, High Intensity Focused Ultrasound, Irreversible Electroporation, and Cryotherapy.

Santiago Uribe-Lewis, MD, PhD

Land Ketsawang, MSc



Contents

Abstract	Page 3
Introduction	Page 4
Focal treatment modalities overview	Page 5
Focal technologies	Page 5
Focal approaches: key differences	Page 8
Methods	Page 9
Literature search strategy	Page 9
Evidence synthesis	Page 9
Results	Page 10
Discussion	Page 13
About the authors	Page 15
References	Page 16

Focaltherapy for localised prostate cancer: systematic review of efficacy and functional outcomes from partial gland ablations with Low Dose Rate Brachytherapy, High Intensity Focused Ultrasound, Irreversible Electroporation, and Cryotherapy.

Abstract

Objective:

To compare oncological and functional outcomes of focal low dose rate I-125brachytherapy (FB),focal high intensity focused ultrasound (HIFU),irreversibleelectroporation (IRE),andfocal cryotherapy for the treatment of localised prostate cancer (PCa).

Methods:

A systematic review of the literature was conducted to obtain clinicalreports on the focal approaches. The searchwas restricted to articles published from May 2015to March 2026. Articles on primary focal monotherapy with peer-review provenance in the Englishlanguage were included. Treatment efficacy outcome was incidence of clinically significant PCa (csPCa; Gleason group 2 or more) recurrence in per protocol control biopsies. Treatment-related effects were assessedfrom patient-reported outcome measures and physician-reported adverse events.

Results:

Ninety-seven articles were appraised including 14 on FB,42 on HIFU, 21on IRE,and 20 on Cryotherapy. Most studies were single-arm prospective analyses. Randomised comparisons between the focal approaches were not identified. In studies with per protocol control biopsies (N=68) the median recurrence of csPCain-field was 3%, 6%, 10%, and 11%after FB, IRE, Cryotherapy and HIFU in 9, 17, 11, and 27 studies respectively. Urinary function reported in 51studies was largely preserved or improved. Sexualfunction assessedin 54 studies was preserved in 35 (65%) studies, and decreased in 12%, 33%, 36% and 55% of FB, HIFU, IRE, and Cryotherapy studies respectively. Bowel function from 18studies was always preserved. The median of any grade 2 adverse event was 5% after Cryotherapy, 8% after HIFU, and 9% after FB and IRE. The median of any grade 3 adverse event was 0% after FB, IRE, and Cryotherapy and 2.8% after HIFU.

Conclusion:

The real-world evidence suggests FBis more efficient at controlling localised PCa relative to HIFU, IRE, and Cryotherapy with a low impact on organ-at-risk function. Randomised comparisons are however needed to confirm long-term efficacy and safety advantages between the focal therapies, whole gland treatment, and active surveillance.

Introduction

Modern prostate magnetic resonance imaging (MRI) provides localisation of tumours together with information on their clinical aggressiveness [1]. This has paved the way for targeted diagnostic approaches and focal treatment where interventions are restricted to areas found to contain neoplasia. By preserving healthy tissues, focal therapy aims to control prostate cancer with a reduced morbidity.

The choice of treatment for localised PCa is a joint decision between the patient and the treating physician. Quality of life after treatment and medical expertise play an important role in treatment selection. However, the risk of tumour progression to more advanced stages takes precedence. For low-risk tumours, patients may be enrolled in an active surveillance program in which close monitoring defers radical treatment, and attendant morbidity, unless the cancer shows signs of progression. Focal therapy may reduce uncertainty stemming from leaving tumours untreated [2]. Although not fully recognised as a standard of care it is increasingly more accepted as an option for intermediate-risk cancer when used in a clinical trial or prospective registry setting (table 1) where patients are closely monitored after treatment.

Focal therapy for prostate cancer has been extensively reviewed in the clinical literature. Valerio et al in 2017 concluded from a systematic synthesis of 37 articles on 7 sources of energy (including brachytherapy, HIFU, IRE, and Cryotherapy) that focal approaches seemed to have a minor impact on genito-urinary function but an undetermined oncological effectiveness due to short follow-up, study heterogeneity, and absence of comparator arms [3]. Ganzer et al in 2018 suggested from an expert consensus that each form of focal therapy may be preferable for different tumour locations, and that tumour risk, size, and prostate volume need consideration [4]. More recent reviews continue to support these initial observations [5-10].

This review captures those clinical studies with an update from the clinical literature to compare oncological efficacy and functional outcomes after FB with those from HIFU, IRE, and Cryotherapy.

Prostate cancer		EAU [12]	AUA/ASTRO [13]	FALCON [14]
Localised	Low-risk	?	⚠	⚠
	Intermediate-risk	?	✓*	✓
	High-risk	✗	✗	✗
Locally Advanced		✗	✗	✗

✓ Recommended ? Trial setting or registry ⚠ AS preferred ✗ Not recommended

*Selected and informed patients, ideally in a prospective study. EAU: European Association of Urology, AUA: American Urology Association, ASTRO: American Society for Radiation Oncology, FALCON: FocaAL Therapy Consensus, AS: Active Surveillance

Table 1. Summary of recommendations for focal therapy in the management of prostate cancer (adapted from [11]).

Focal treatment modalities overview

Low dose rate brachytherapy

Low dose rate (LDR) brachytherapy seeds emit ionising radiation that breaks molecular bonds in DNA resulting in cellular death by apoptosis. Each seed is a titanium tube 4.5 mm long by 0.8 mm in diameter (similar to a grain of rice) that encapsulates the radioisotope, most commonly Iodine-125 (I-125, half-life of 59.6 days) or Palladium-103 (Pd-103, half-life of 17 days). The number and position of seeds is planned using dedicated software to achieve the prescription dose of 145 Gy for I-125 to the target area identified from MRI and biopsy information. The seeds are permanently implanted via the perineum under real-time transrectal ultrasound (TRUS) guidance. After the implant, seed placement is verified by computerised tomography or MRI which enable a record of the doses delivered to the target area(s) and surrounding organs. See figure 1 and table 2 for further technical detail and comparisons on the selected focal modalities.

High-intensity focused ultrasound (HIFU)


Ultrasound frequency of 0.8–7 MHz transforms into mechanical energy generating temperatures between 65°C and 90°C. The thermal ablation causes protein denaturation and coagulative necrosis. Once dimensions of the prostate are acquired and the target area identified, a 'cigar shaped' ablation volume (1.5–2.0 mm wide by 15–20 mm long covering 50 to 300 mm³) is guided to the target area in real-time by means of dedicated software. Fusion of mpMRI with ultrasound images assists with lesion location, and robotic technology enhances precision. The procedure is non-invasive as energy is delivered from the ultrasound probe housed in a cooling balloon within the rectum. In large prostate glands delivery of energy to the anterior prostate may not be feasible due to the fixed focal length.

Irreversible electroporation (IRE)

Electroporation forms nanopores in the cellular membrane that become permanent leading to loss of homeostatic regulation and cell death. Paired electrodes of anode (P+) and cathode (P-) needles are placed transperineally under TRUS guidance to deliver a series of short-duration high-voltage electrical pulses. The electrodes are placed 10 mm to 20 mm apart to create an ablation zone that surrounds the lesion. Multiple pairs of electrodes can be sequenced to increase the size of the ablation zone (see figure 2). The active tip (exposed length) of the electrodes determines the extent of the area to be treated.


Cryotherapy

Cryotherapy involves two cycles of freezing and thawing. Lethal freezing temperatures (below -40°C) damage the cellular membrane and microvascular endothelium causing coagulative necrosis. Thawing creates a hypotonic environment that leads to cell swelling, further membrane damage, and oedema. Thin cryotherapy needles (cryoprobes) are inserted via the perineum. Compressed Argon gas expands within the tip of the cryo-needle to form an ice ball and helium gas enables the thawing phase. The length and diameter of the cryoprobes determine the shape and size of the ice ball which is monitored under TRUS guidance. Thermocouples are used to monitor temperatures in organs at risk and a warmed (40°C) balloon is placed in the rectum to control rectal temperature.



a. FB:


Low dose rate brachytherapy seeds permanently implanted in and around the target area. The image shows implants of seeds in strands. A treatment plan based on the size and location of the tumour, and proximity to organs at risk, determines strand configurations.



b. HIFU:


The probe emits ultrasound from the rectum with a fixed focal length, hence requires movement and rotation to reach targets in different areas of the prostate, and compresses the prostate to reach tumours in the anterior part of the gland (prostate deformation shown in blue dashed line).

Figure 1 Page 6-7. Illustrations of single lesion focal treatment, example shown with lesion at mid anterior and right posterior base for each modality to highlight possible treatment approach at different locations around key organs; a. Focal Brachytherapy (FB), b. High Intensity Focused Ultrasound (HIFU), c. Irreversible electroporation (IRE), d. Cryotherapy. Shown in transverse and sagittal orientation with estimated treatment / energy margins.



c. IRE:

Four and six electrodes implanted into the prostate. Each probe is paired and sequenced e.g. a and d, b and e, c and f, etc., to generate a high voltage field around the tumour. The exposed length of the electrodes are used to match to the length and shape of the tumour.



d. Cryotherapy:

Thin needles (cryoprobes) expand Argon gas at the tip to create an 'ice ball' in the target area. The cryoprobe tip size determines the volume/shape and isotherm dimensions. Isotherm dimensions vary on the needle used, with temperatures reaching -400°C , -200°C , and 0°C , within the white, light blue, and dark blue dashed lines, respectively.

	FB	HIFU	IRE	Cryotherapy
Treatment principle	Apoptosis by ionising radiation of DNA ¹	Apoptosis by heat/ concentrated mechanical frequency ²	Permanent nanopores lead to loss of homeostasis, may trigger immune response ³	Apoptosis and necrosis by tissue freeze/thaw ⁴
Energy delivery	145Gy, 59.6 days half-life ⁵	200kHz and 7MHz generating 60°C to 90°C or Hyperchoic 20 to 48 Watts ²	1500 V/cm at 90-microsecond monopolar pulses ⁵	38°C to -40°C time-controlled temperature gradient ⁴
Selection criteria ⁱ	<ul style="list-style-type: none"> • Good urinary function (high IPSS)⁷ • Prostate <20cc OR >60cc (100cc is technically feasible) with minimal pubic interference¹¹ • Prior TURP possible⁷ 	Prostate >40cc due to limited focal length ⁸	<ul style="list-style-type: none"> • Pacemakers must be identified⁹ • Prostate volume 65cc with minimal pubic arch interference¹² 	<ul style="list-style-type: none"> • Prostate size between 40-60cc¹⁰ • Men with prior TURP treatment not suitable¹³
Patient preparation ⁱⁱ	Separate TRUS volume study needed for pre-loaded kit ordering	<ul style="list-style-type: none"> • Dorsal lithotomy or right lateral decubitus • TURP may be performed to reduce prostate volume and minimise necrotic debris from procedure¹⁵ • Immobilisation critical¹⁷ • Rectal cooling from balloon housing the probe¹⁷ • Removal of gas bubbles desired to ensure uniform wave propagation² 	Muscle relaxant to prevent contractions ¹⁶	<ul style="list-style-type: none"> • Warmed urethral catheter, non-target tissue thermal protection¹⁸ • Circulating heated rectal balloon¹⁹ • Thermocouples placed at OAR²⁰
Imaging and planning ⁱⁱⁱ	<ul style="list-style-type: none"> • TRUS real-time • Intra-operative real-time analysis after seed each implant²⁴ 	<ul style="list-style-type: none"> • mpMRI fusion with TRUS real-time²¹ • Ultrasound and HIFU probe are combined as one 	mpMRI fusion with TRUS real-time ²²	mpMRI fusion with TRUS real-time ²³
Targeting & instrument placement ^{iv}	<ul style="list-style-type: none"> • Seeds placed in and around lesion²⁴ • 13 brachytherapy needles (median) inserted transperineally under TRUS²⁴ • 39 permanent implant seeds with the option of being encased/linked and spaced (median)²⁴ 	<ul style="list-style-type: none"> • Proximal energising of lesion via software guidance² • Non-invasive • 50-56 mm dia visualization and treatment probe • Robotic arm to direct treatment probe within the rectum or mounted on articulated arm and fixed into position 	<ul style="list-style-type: none"> • Lesion must be between paired probes²⁵ • 4-6 probes parallel configuration inserted transperineally under TRUS²⁵ • 10-20 mm spacing⁶ 	<ul style="list-style-type: none"> • Ice-ball formed within and around the lesion • 1-4 cryotherapy needles placed within 20 mm of each other inserted transperineally under TRUS²⁶ • Needle dimensions reflect ice-ball shape • Temperature probes in urethra sphincter, rectal wall, NVB²⁷
Treatment margin	7 mm isotropic expansion of the focal-gross tumour volume (F-GTV) ⁵	<ul style="list-style-type: none"> • 10 mm beyond visible lesion²⁸ • Margin beyond the prostate capsule when disease abuts the capsule are recommended³⁰ 	5-7 mm beyond the visible lesion ²⁹	A minimum 5 mm of <20°C around the target lesion/ zone ²⁷
Margin: Safety from critical structures ^v	<ul style="list-style-type: none"> • 3 mm from urethra²⁴ • Urethra volume 0.1cc (V200%) receiving max dose of 164.6 Gy (mean)⁵ • Rectum volume 0.05cc (V100%) receiving max dose 95.8 Gy (mean)⁵ • Peri-rectal spacer placed after implant absorbed after 6 months reduces rectal dose⁵ 	5 mm from rectum and urinary sphincter ¹⁷	<ul style="list-style-type: none"> • 5 mm from critical structure⁶ • 3 mm from prostate capsule⁶ 	<ul style="list-style-type: none"> • Urethra temperature to be maintained at 38-40°C¹⁸ • Rectum target temp 40°C with warmer cryotherapy needles placed between the prostate and rectum, and/or rectal lumen washed with 40°C water for -10°C rectal temps¹⁹ • NVB not below -35°C³¹
Procedure time	36 minutes (mean), real-time implant ²⁴	63 minutes, treatment area of 14.67mL (hemi-ablation) in 7 slices ¹⁷	60 minutes ¹⁶	70 minutes (median), freezing thawing 6-7 minutes in 2 cycles ²⁷
Post-procedure patient care activity ^{vi}	<ul style="list-style-type: none"> • Cystoscopy post implant, visual seed check within urethra or bladder⁵ • Urethral catheterisation for urinary retention, occurs in <2% of whole gland therapy procedures⁵ 	Bladder catheter for 5-7 days ¹⁷	Urine catheter in place for 2 – 10 days depending on the size of the ablation zone and proximity to intraprostatic urethra ¹⁶	<ul style="list-style-type: none"> • Urethral catheter 3-7 days³² • Belladonna / opium suppositories²⁷
Treatment quality assurance	CT or mpMRI post-implant dosimetry ⁵	Contrast enhanced TRUS correlates with mpMRI and biopsy histology ³³	Contrast enhanced ultrasound for perfusion defect in relation to ablation zone ³⁴	Visual confirmation of treated area with a follow up scan.

Common across modalities:

- Patient exclusion include rectal fistulae and stenosis, calcification, history of abdominal perineal resection for rectal cancer, or other major rectal pathology including inflammatory bowel disease. Lesion sitting directly at the posterior apex of the prostate, near the external urinary sphincter, face restrictions. Risk of permanent urinary incontinence increased if ablation zone significantly overlaps with the sphincter.
 - Patient preparation includes antibiotic prophylaxis, dorsal lithotomy positioning, general anaesthesia, catheter placement (or aerosolised jelly) for landmarking and drainage, complete patient mobilization is non-critical
 - All tumours (index lesions) and critical structures mapped with software enabling targeted treatment and real-time ultrasound monitoring
 - Procedure is minimally invasive with template and stepper used for alignment under real-time guidance
 - Margins at the discretion of the physicians.
 - Alpha-blocker medication ranges from 6 months depending on modality
- FB: Focal Brachytherapy; HIFU: High-intensity Focal Ultrasound; IRE: Irreversible Electroporation; Gy: Gray; TURP: Transurethral Resection of the Prostate; OAR: Organs At Risk; TRUS: Transrectal Ultrasound; mpMRI: multi-parametric Magnetic Resonance Imaging; NVB: Neurovascular Bundle; F-GTV: focal-gross tumour volume; CT: Computed Tomography.

See references on page 19.

Table 2. Summary of focal approaches highlighting key differences

Methods

Literature search strategy

The PubMed and EMBASE databases were queried for literature on focal approaches for the treatment of prostate cancer with the following PICO terms:

Population – Prostate cancer

Intervention – Focal therapy

Comparator – Energy sources

Outcome – Efficacy (oncological control) and Functional (treatment-related effects).

To obtain a broad view an initial search explored availability of reviews or clinical trials on focal therapy for prostate cancer within the last 10 years (figure 2) from 25 May 2025. The reviews guided hand searches for selection of relevant clinical studies. A second search was conducted to obtain up to date clinical studies from 25 May 2025 to 10 March 2026.

Evidence synthesis

Data from each study was collected in Excel and summarized within R statistical environment (version 4.5.1) [15] using readxl [16], tidyverse [17], gtsummary [18], and RColorBrewer [19] packages. Age means, when not reported, were obtained from available summary statistics using the Box-Cox method in 'estmeansd' [20].

Efficacy outcomes, in-field and out-field recurrence, were defined as the proportion of any Gleason group detected in per protocol control biopsies. Clinically significant recurrence when Gleason group 2 or more was identified. Functional outcome impact was defined by evidence of a statistically significant change from baseline of PROMs (IPSS, AUS-SS, EPIC-26, IIEF-5) beyond an acute phase (>3 months after treatment).




Figure 2. Overview of literature search strategy. The first search was for reviews in the last 10 years up to May 2025, and a second search for studies from May 2025 to March 2026.

Results

The literature search found 97 studies pertaining to focal therapy of localised prostate cancer with the selected energy sources (figure 2 and table 3). For FB 14 studies encompassed 865 cases with a median of 30 cases per study [21-34]. The largest body of evidence was for HIFU with 42 studies on 5,246 cases and a median of 73 cases per study [35-76]. For IRE 21 studies included 1,763 cases and a median of 50 cases per study [77-97]. For Cryotherapy there were 20 studies on 2,424 cases and a median 76 cases per study [98-116,71]. One study reported on HIFU and Cryotherapy [71].

Most reports were from prospective registries (~70%) or retrospective analyses (~30%) (table 3). The single randomised trial found compared focal versus extended IRE without a significant difference in outcomes [96]. Patient age was similar between all treatment modalities, and most studies treated men with prostate cancer of a low and intermediate risk

of relapse. Inclusion of men with high-risk disease was more prevalent in reports on cryotherapy (55% of studies) and HIFU (33.8% of studies). The length of patient follow-up was short overall with FB studies reporting the longest follow up length at a median of 36 months (table 3). For all therapies the number of patients per study grew over time (figure 3).

Clinically significant PCa recurrence detected in treated (in-field) and untreated (out-field) areas in per protocol biopsies was reported by 64 and 61 studies respectively (figure 3 and table 4). The median in-field recurrence rate was 3% after FB (from 233 biopsies across 9 studies), 11% after HIFU (from 1,953 biopsies in 27 studies), 6% after IRE (from 1,108 biopsies in 17 studies), and 10% after Cryotherapy (from 761 biopsies across 11 studies). The median out-field recurrence was 3% for FB (in 9 studies), 8% for HIFU (in 27 studies), 9% after IRE (in 16 studies), and 12% after Cryotherapy (across 9 studies).

	FB	HIFU	IRE	Cryo
Study characteristic	N = 14 ¹	N = 42 ¹	N = 21 ¹	N = 20 ¹
Number of cases	30 (24, 51)[865]	73 (44, 113)[5,246]	50 (25, 106)[1,763]	76 (53, 127)[2,424]
Study type				
Prospective	10 (71%)	33 (79%)	15 (71%)	13 (65%)
Retrospective	4 (29%)	9 (21%)	5 (24%)	7 (35%)
RCT	0 (0%)	0 (0%)	1 (4.8%)	0 (0%)
Focal approach				
Ultrafocal	1 (7.1%)	0 (0%)	0 (0%)	0 (0%)
Focal	8 (57%)	16 (38%)	20 (95%)	11 (55%)
Focal or Partial	1 (7.1%)	13 (31%)	1 (4.8%)	1 (5.0%)
Partial or Hemigland	4 (29%)	13 (31%)	0 (0%)	8 (40%)
Age means (years)	67 (63, 71)	66 (65, 68)	66 (64, 67)	67 (65, 69)
Unknown	0	2	0	1
Risk groups				
Low	0 (0%)	2 (4.8%)	0 (0%)	1 (5.0%)
Int	1 (7.1%)	3 (7.1%)	1 (4.8%)	4 (20%)
Low-Int	12 (86%)	23 (55%)	16 (76%)	4 (20%)
Low-Int-High	1 (7.1%)	12 (29%)	4 (19%)	6 (30%)
Int-High	0 (0%)	2 (4.8%)	0 (0%)	5 (25%)
Followup length median (months)	36 (24, 65)	20 (12, 30)	19 (12, 36)	27 (15, 36)
Control biopsy reason				
Per protocol	8 (57%)	30 (71%)	19 (90%)	12 (60%)
For cause	3 (21%)	9 (21%)	2 (9.5%)	7 (35%)
None	3 (21%)	3 (7.1%)	0 (0%)	1 (5.0%)

¹Median (Q1,Q3) [Sum]; n (%); Median (Q1,Q3)

Table 3.





Figure 3. Overall data structure and recurrence. Top panel shows all 97 studies by year of publication and patient numbers with circle sizes by follow-up lengths. Lines represent slopes from linear regression models. Bottom panels show clinically significant (csPCa) in-field and out-field recurrence in studies with per protocol biopsies.

	FB	HIFU	IRE	Cryo
Perprotocol biopsies	N = 9 ¹	N = 30 ¹	N = 19 ¹	N = 12 ¹
Number of biopsied men	26 (17,34) [233]	57 (31,90) [1,953]	40 (21, 101)[1,108]	51(39, 76) [761]
In-field positive - any Gleason (%) ²	4 (0, 10) [0:17]	15(10, 23) [0:40]	12(6, 18)[0:37]	11(7, 20) [2:37]
In-field positive - csPCa(%) ³	3 (0, 7) [0:10]	11(5, 15) [0:23]	6 (4, 9) [0:14]	10 (6, 13)[2:23]
Unknown	0	3	2	1
Out-field positive - any Gleason (%) ²	12(7, 18) [0:41]	16(9, 24) [4:38]	16(6, 27) [0:50]	17(11,23) [1:40]
Unknown	0	0	0	2
Out-field positive - csPCa(%) ³	3 (0, 8) [0:16]	8 (5, 10) [0:38]	9 (6, 15) [0:17]	12(9, 13) [1:22]
Unknown	0	3	3	3

1Median (Q1,Q3)[Sum]; Median (Q1,Q3)[Min:Max]
 2Inmen with per protocol biopsy (includesbiopsieswith evidenceof treatment effect after brachytherapy).
 3Gleasongroup 2 or more.

Table 4. Summary of oncological outcomes.

Most studies that used PROMs for urinary function showed no impact on urinary symptoms (figure 4). Two studies on HIFU [35, 76], 2 on IRE[85, 92], and 3 on cryotherapy [102, 106, 109] reported better urinary symptoms post-ablation. One study showed worse urinary function after FB (6 months after implantation)[25]. The impact on patient-reported sexual function was more pronounced across all treatments; symptoms were worse in 1 FB study (12%, and clinically significant only in men with erectile disfunction at baseline at 6 months and 3 years)[21], 7 HIFU studies (33%)[46, 48, 54, 56, 57, 73, 74], 5 IRE studies (36%)[77, 78, 85, 88, 93], and 6 cryotherapy studies (55%) [100, 102, 104, 106, 108, 109]. Improvement in sexual function was never reported. An impact on bowel function was assessed in 18 studies all of which showed no significant change in patient-reported bowel symptoms.

Grade 2 and 3 adverse events were reported in 54 and 52 studies respectively (figure 5). Genitourinary (GU) events were in general low across therapies; for grade 2 GU events the median was 5% in FB, 6% after Cryotherapy, and 9% for HIFU and IRE with the top quartile below 20%. A notable outlier was one study on HIFU with 62% (18 events in 29 patients) [42]. For grade 3 GU events the median was 0% after FB, IRE, and Cryotherapy, and 4% after HIFU with the top quartile below 10%. Gastrointestinal (GI) events were rare with most grade 2 and 3 events below 1% and 0.5% of treated men respectively. A notable outlier was one study of FB with grade 2 GI of 6.6% (2 events in 30 patients)[25].




Figure 4. Focal therapy-related toxicity. Urinary, sexual, and bowel function was assessed by 51, 54, and 18 studies, respectively, that compared PROMs before and after treatment beyond a 3 months acute phase. Shown above bars are the percentage of studies with an impact on function (impact defined as a statistically significant difference relative to pre-treatment values). The number of studies are indicated inside the bars.

Discussion

The field of focal therapy continues to gain momentum with its place as a treatment option for localised prostate cancer increasingly more established [12-14]. A growing body of evidence is beginning to support the promise of equal efficacy to whole-gland radical treatment with reduced impact on nearby organ function. Hence men may find an offer of focal therapy more appealing than leaving tumours untreated, whilst accepting a risk of leaving healthy-looking areas untreated, if their quality of life is preserved. Focal therapy has also shown a favourable therapeutic ratio (effective cancer control with low morbidity) as salvage treatment [117].

Level 1 evidence is needed to compare outcomes between the different focal therapy approaches, whole gland treatment, and active surveillance. The

outcomes herewith reported should be interpreted with caution. Most studies found were from single-arm prospective registries with relatively small case numbers. Patient selection criteria were not universal. Most studies assessed outcomes from mixed low and intermediate-risk cases, and several studies included men with a high risk of relapse in varying proportions; notably more than half of IRE studies included 13% to 34% high-risk cases (not shown). The ablation strategies varied from ultrafocal to partial or hemigland and several studies reported outcomes from mixed populations of focal and partial (up to 31% of HIFU studies). This has implications for both efficacy and functional outcomes. Interpretation of functional outcomes is further hampered by variation in the type of patient-reported instruments used, follow up protocols, and applied statistics.




Figure 5. Adverse events. Genitourinary (GU) and gastrointestinal (GI) events are presented as the proportion of patients with an event in the study population. The number of studies reporting adverse events for each focal therapy approach is shown inside parenthesis.

This review found 97 studies on 4 sources of energy for primary focal therapy (I-125 Brachytherapy, HIFU, IRE, and Cryotherapy). A recent systematic review by Hopstaken et al published in 2022 found 72 studies on 8 sources of energy. At that time rates of in-field clinically significant cancer were reported as 0% for FB (LDR or HDR), 14.7% for HIFU, 8.5% for IRE, and 15% for cryotherapy [5]. Here, with a larger body of evidence, we obtained a similar pattern of efficacy outcomes between the sources of energy. The length of follow up also remains short which continues to challenge an assessment of the oncological effectiveness of focal therapies.

Seventy studies conducted per protocol biopsies. From these, I-125 FB showed the lowest rates of in-field and out-field recurrence. Recurrence after brachytherapy would be further reduced if presence of a radiation treatment effect in the biopsy specimen was not considered a recurrence. The strength of evidence, in terms of number of studies and study participants, was however lower for FB studies. Lower proportions of high-risk participants in the FB studies may account for the lower recurrence rates observed. Of note, exclusion of studies with high-risk men continued to show a similar pattern with median recurrence rates of 3%, 4%, 7%, and 12% for FB, cryotherapy, IRE, and HIFU, respectively (not shown).

Per protocol post-treatment control biopsies were taken from 6 months to 3 years post-ablation (in 26 studies within the first year and in 43 studies between year 1 and 3). These short-term measures of efficacy may need to be reconsidered once longer-term outcomes become available. Indeed, a number of studies abandoned 6-month biopsies due to very low incidences of significant disease. The FOCAL consensus has recommended control biopsies should be performed within the first 12 months post-treatment [14].

The impact of focal therapy on urinary, sexual, and bowel function was assessed by 51, 54, and 18 studies, respectively, that compared PROMs (IPSS, AUA-SS, IIEF-5, EPIC) before and after treatment (beyond an acute phase of 3 months). Urinary function was largely preserved, and in some instances improved. Only one study reported worse urinary function (increased IPSS) 6 months after focal/partial brachytherapy (N=30 patients), although it was significantly better than that seen in whole-gland implant controls (N=30 patients) [25]. Sexual function was more frequently negatively impacted but for the most part preserved; 65%

(35/54) of studies reported no difference in sexual function (IIEF-5 and/or EPIC) relative to baseline. No study showed improvement in sexual function. Bowel function, although much less frequently assessed, was always preserved. Rectal spacers were used in one study [21].

For this review the impact of organ-at-risk function was defined by statistical significance. This may not necessarily translate into a clinically meaningful difference. The type of instrument used to quantify patient-reported symptoms may have a considerable bearing on outcome interpretation. For example, in three studies that reported a decline in sexual function there was discrepancy on statistical significance between SHIM and EPIC-SF scores at 3 years in 39 patients after HIFU [54], at 3 years in 50 patients after HIFU [74], and at 1 year on 50 patients after Cryotherapy [109]. The latter study showed the baseline score was the only significant predictor of sexual function scores after treatment. Larger study populations and longer-term follow up may be necessary to resolve these discrepancies.

Adverse event classification systems included Clavien-Dindo (38 studies), CTCAE (Common Terminology Criteria for Adverse Events) (14 studies), and RTOG (Radiation Therapy Oncology Group) morbidity scale (1 study). Adverse events were low across the 4 therapies. Two notable outliers were Duwe et al with 62% grade 2 GU events after HIFU [42] and Kim et al with 6.6% grade 2 GI events after FB [25]. The former used Clavien-Dindo and the latter the only study that used RTOG criteria for rectal toxicity. In both systems grade 2 indicates requirement for pharmacological treatment. Grade 3 GU and GI events were more frequent after HIFU but remained below 15% and 1% respectively. The timing of adverse event occurrence was only reported by one third of the studies.

In conclusion, the outcomes reported herewith on efficacy and function support findings from recent literature reviews and meta-analyses collectively comparing focal therapy approaches [5-11]. The real-world evidence suggests FB is more efficient at controlling localised PC relative to HIFU, IRE, and Cryotherapy with a low impact on organ-at-risk function. However, considerable heterogeneity in study protocols preclude robust conclusions. Randomised comparisons between focal therapies, whole gland treatment, and active surveillance are needed to confirm the emerging consensus on a favourable therapeutic ratio of partial gland ablation.

About the authors

Santiago Uribe-Lewis is director of BXTConsultants Ltd since 2017. He obtained his degree in Medicine from the Universidad del Rosario (Bogotá, Colombia) in 1995, and PhD in Molecular Biology/Epigenetics from Imperial College (London, UK) in 2006, followed by post-doctoral research on Cancer Epigenetics at the University of Cambridge (UK). In 2014 he joined BXTConsultants Ltd and from 2016 an honorary appointment with the Department of Urology at the Royal Surrey NHS Foundation Trust (Guildford, UK).

Land Ketsawang is International Marketing Director of BXTA since 2023. He obtained his Bachelor Engineering degree in Computer Aided Engineering Design from Sheffield Hallam University (UK) in 2002, then completed a Masters in Biomedical Engineering & Medical Physics from University College London (UK) in 2006. He has worked in the medical technology industry since 2009 with engineering roles at Siemens-Healthineers MRI Technology (Oxford, UK), then moving into EMEA commercial roles at BD Life Sciences (UK and Switzerland), before joining BXTA (UK).

References

1. Turkbey B, Rosenkrantz AB, Haider MA et al. Prostate Imaging Reporting and Data System Version 2.1:2019 Update of Prostate Imaging Reporting and Data System Version 2. *Eur Urol* 2019; 76: 340-51.
2. Taylor KL, Luta G, Hoffman RM et al. Quality of life among men with low-risk prostate cancer during the first year following diagnosis: the PREPARE prospective cohort study. *Transl Behav Med* 2018; 8: 156-65.
3. Valerio M, Cerantola Y, Eggener SE et al. New and Established Technology in Focal Ablation of the Prostate: A Systematic Review. *Eur Urol* 2017; 71: 17-34.
4. Ganzer R, Arthanareeswaran VKA, Ahmed HU et al. Which technology to select for primary focal treatment of prostate cancer? - European Section of Urotechnology (ESUT) position statement. *Prostate Cancer and Prostatic Diseases* 2018; 21: 175-86.
5. Hopstaken JS, Bomers JGR, Sedelaar MJP et al. An Updated Systematic Review on Focal Therapy in Localized Prostate Cancer: What Has Changed over the Past 5 Years? *Eur Urol* 2022; 81: 5-33.
6. Busby D, Rich JM, Grauer R et al. Bimodal and Erectile Functional Outcomes of Partial Prostate Ablation: A Systematic Review and Meta-analysis of Prospective Studies. *Urology* 2023; 182: 14-26.
7. Tay KJ, Fong KY, Stabile A et al. Established focal therapy-HIFU, IRE, or cryotherapy-where are we now?-a systematic review and meta-analysis. *Prostate Cancer Prostatic Dis* 2024.
8. Mohamad O, Nicosia L, Mathier E et al. Focal brachytherapy as definitive treatment for localized prostate cancer: A systematic review and meta-analysis. *Brachytherapy* 2024; 23: 309-20.
9. Piramide F, Veccia A, Tzelvels L et al. Sexual function outcomes in men undergoing minimal invasive ablative techniques for prostate cancer: a ESRU/YAUrotech systematic review and pooled analysis. *Minerva Urol Nephrol* 2025; 77: 285-97.
10. Le Guevelou J, Peyrottes A, Barret E et al. Focal radiation therapies for localized prostate cancer: a systematic review of clinical evidence and outcomes. *Prostate Cancer Prostatic Dis* 2026; 29: 303-11.
11. Roldan-Testillano R, Rodriguez-Sanchez L, Covarrubias C et al. Advances in focal therapy for prostate cancer: current modalities, outcomes, and future directions. *Prostate Int* 2026; 14: 1-9.
12. EAU Guidelines. <https://uroweb.org/guidelines/prostate-cancer>.
13. Eastham JA, Barocas DA, Chu CE et al. Clinically Localized Prostate Cancer: AUA/ASTRO Guideline Amendment (2026). *J Urol* 2026; 101097JU00000000000005060.
14. Rodriguez-Sanchez L, Reiter R, Rodriguez A et al. The FocAL therapy CONsensus (FALCON): enhancing partial gland ablation for localized prostate cancer. *BJU Int* 2024; 134: 50-3.
15. R Core Team (2025). *R: A Language and Environment for Statistical Computing*. R Foundation for Statistical Computing, Vienna, Austria. <https://www.R-project.org/>
16. Wickham H, Bryan J (2025). *readxl: Read Excel Files*. doi:10.32614/CRAN.package.readxl <<https://doi.org/10.32614/CRAN.package.readxl>>, R package version 1.4.5, <https://CRAN.R-project.org/package=readxl>
17. Wickham H, Averick M, Bryan J et al. Welcome to the Tidyverse. *Journal of Open Source Software* 2019; 4.
18. Sjoberg DD, Whiting K, Curry M, Lavery JA, Larmarange J. Reproducible Summary Tables with the gtsmmary Package. *The R Journal* 2021; 13.
19. Neuwirth E (2022). *RColorBrewer: ColorBrewer Palettes*. doi:10.32614/CRAN.package.RColorBrewer <https://doi.org/10.32614/CRAN.package.RColorBrewer>, R package version 1.1-3, <https://CRAN.R-project.org/package=RColorBrewer>.
20. McGrath S, Zhao X, Steeler R, Benedetti A (2023). *estmeansd: Estimating the Sample Mean and Standard Deviation from Commonly Reported Quantiles in Meta-Analysis*. doi:10.32614/CRAN.package.estmeansd <https://doi.org/10.32614/CRAN.package.estmeansd>, R package version 1.0.1, <https://CRAN.R-project.org/package=estmeansd>.
21. Adhami M, Cheng J, Anderson E et al. Functional outcomes and toxicity following focal low-dose-rate brachytherapy for low-intermediate risk prostate cancer. *BJUI Compass* 2026; 7: e70119.
22. Adhami M, Cheng J, Anderson E et al. Oncological outcomes post focal low-dose-rate brachytherapy in low-intermediate risk prostate cancer. *BJUI Compass* 2026; 7: e70129.
23. Anderson E, Smyth LML, O'Sullivan R et al. Focal low dose-rate brachytherapy for low to intermediate risk prostate cancer: preliminary experience at an Australian institution. *Transl Androl Urol* 2021; 10: 3591-603.
24. Graff P, Portalez D, Lusque A et al. IDEAL2a Phase II Study of Ultrafocal Brachytherapy for Low- and Intermediate-risk Prostate Cancer. *Int J Radiat Oncol Biol Phys* 2018; 102: 903-11.
25. Kim TH, Kim JN, Yu YD et al. Feasibility and early toxicity of focal or partial brachytherapy in prostate cancer patients. *J Contemp Brachytherapy* 2020; 12: 420-6.
26. King MT, Nguyen PL, Boldbaatar N et al. Long-term outcomes of partial prostate treatment with magnetic resonance imaging-guided brachytherapy for patients with favorable-risk prostate cancer. *Cancer* 2018; 124: 3528-35.
27. Kunogi H, Wakumoto Y, Kawamoto T et al. Focal low-dose-rate prostate brachytherapy for low- and intermediate-risk prostate cancer. *J Contemp Brachytherapy* 2020; 12: 554-61.
28. Langley S, Uribe J, Uribe-Lewis S et al. Hemi-ablative low-dose-rate prostate brachytherapy for unilateral localised prostate cancer. *BJU Int* 2020; 125: 383-90.
29. Langley S, Uribe J, Uribe-Lewis S et al. Is hemi-gland focal LDR brachytherapy as effective as whole-gland treatment for unilateral prostate cancer? *Brachytherapy* 2022; 21: 870-6.
30. Mahdavi SS, Spadinger IT, Salcudean SE et al. Focal application of low-dose-rate brachytherapy for prostate cancer: a pilot study. *J Contemp Brachytherapy* 2017; 9: 197-208.
31. Matsuoka Y, Uehara S, Toda K et al. Focal brachytherapy for localized prostate cancer: 5.7-year clinical outcomes and a pair-matched study with radical prostatectomy. *Urol Oncol* 2022; 40: 161e15-e23.
32. Saito K, Matsuoka Y, Toda K et al. Medium-term oncological and functional outcomes of hemi-gland brachytherapy using iodine-125 seeds for intermediate-risk unilateral prostate cancer. *Brachytherapy* 2021; 20: 842-8.
33. Srougi V, Barret E, Nunes-Silva I et al. Focal brachytherapy for localized prostate cancer: Urinary toxicity depends on tumor location. *Brachytherapy* 2017; 16: 988-92.
34. Ta MH, Nunes-Silva I, Barret E et al. Focal Brachytherapy for Localized Prostate Cancer: Midterm Outcomes. *Pract Radiat Oncol* 2021; 11: e477-e85.
35. Abreu AL, Peretsman S, Iwata A et al. High Intensity Focused Ultrasound Hemigland Ablation for Prostate Cancer: Initial Outcomes of a United States Series. *J Urol* 2020; 204: 741-7.
36. Ahmed HU, Dickinson L, Charman S et al. Focal Ablation Targeted to the Index Lesion in Multifocal Localised Prostate Cancer: a Prospective Development Study. *Eur Urol* 2015; 68: 927-36.
37. Albinini S, Aoun F, Bellucci S et al. Comparing High-Intensity Focal Ultrasound Hemiablation to Robotic Radical Prostatectomy in the Management of Unilateral Prostate Cancer: A Matched-Pair Analysis. *J Endourol* 2017; 31: 14-9.
38. Annot A, Olivier J, Valtille P et al. Extra-target low-risk prostate cancer: implications for focal high-intensity focused ultrasound of clinically significant prostate cancer. *World J Urol* 2019; 37: 261-8.
39. Bacchetta F, Martins M, Regusci S et al. The utility of intraoperative contrast-enhanced ultrasound in detecting residual disease after focal HIFU for localized prostate cancer. *Urol Oncol* 2020; 38: 846 e1-e7.

40. Bass R, Fleschner N, Finelli A et al. Oncologic and Functional Outcomes of Partial Gland Ablation with High Intensity Focused Ultrasound for Localized Prostate Cancer. *J Urol* 2019; 201: 113-9.
41. Dellabella M, Branchi A, Di Rosa M et al. Oncological and functional outcome after partial prostate HIFU ablation with Focal-One(R): a prospective single-center study. *Prostate Cancer Prostatic Dis* 2021; 24: 1189-97.
42. Duwe G, Boehm K, Haack M et al. Single-center, prospective phase 2 trial of high-intensity focused ultrasound (HIFU) in patients with unilateral localized prostate cancer: good functional results but oncologically not as safe as expected. *World J Urol* 2023; 41: 1293-9.
43. Ehdai B, Tempany CM, Holland F et al. MRI-guided focused ultrasound focal therapy for patients with intermediate-risk prostate cancer: a phase 2b, multicentre study. *Lancet Oncol* 2022; 23: 910-8.
44. Ehrengut C, Alberalar ND, Franz T et al. Evaluation of the TARGET score for MRI-based detection of prostate cancer recurrence in patients treated with HIFU hemiablation. *Eur J Radiol* 2026; 195: 112572.
45. Feijoo ER, Sivaraman A, Barret E et al. Focal High-intensity Focused Ultrasound Targeted Hemiablation for Unilateral Prostate Cancer: A Prospective Evaluation of Oncologic and Functional Outcomes. *Eur Urol* 2016; 69: 214-20.
46. Ganzer R, Hadaschik B, Pahernik S et al. Prospective Multicenter Phase II Study on Focal Therapy (Hemiblation) of the Prostate with High Intensity Focused Ultrasound. *J Urol* 2018; 199: 983-9.
47. Ghai S, Finelli A, Corr K et al. MRI-guided Focused Ultrasound Focal Therapy for Intermediate-Risk Prostate Cancer: Final Results from a 2-year Phase II Clinical Trial. *Radiology* 2024; 310: e231473.
48. Glybochko PV, Amosov AV, Krupinov GE, Petrovskii NV, Lumpov IS. Hemiblation of Localized Prostate Cancer by High-Intensity Focused Ultrasound: A Series of 35 Cases. *Oncology* 2019; 97: 44-8.
49. Guillaumier S, Peters M, Arya M et al. A Multicentre Study of 5-year Outcomes Following Focal Therapy in Treating Clinically Significant Nonmetastatic Prostate Cancer. *Eur Urol* 2018; 74: 422-9.
50. Hamdy FC, Elliott D, le Conte S et al. Partial ablation versus radical prostatectomy in intermediate-risk prostate cancer: the PART feasibility RCT. *Health Technol Assess* 2018; 22: 1-96.
51. Hoquetis L, Malavaud B, Game X et al. MRI evaluation following partial HIFU therapy for localized prostate cancer: A single-center study. *Prog Urol* 2016; 26: 517-23.
52. Johnston MJ, Emara A, Noureldin M, Bott S, Hindley RG. Focal High-intensity Focused Ultrasound Partial Gland Ablation for the Treatment of Localised Prostate Cancer: A Report of Medium-term Outcomes From a Single-center in the United Kingdom. *Urology* 2019; 133: 175-81.
53. Kaneko M, Ramacciotti LS, Inoue Y et al. Impact of magnetic resonance imaging visibility of prostate cancer on partial gland ablation. *BJUI Compass* 2025; 6: e70065.
54. Kaufmann B, Raess E, Schmid FA et al. Focal therapy with high-intensity focused ultrasound for prostate cancer: 3-year outcomes from a prospective trial. *BJU Int* 2024; 133: 413-24.
55. Ladjevardi S, Ebner A, Femic A et al. Focal high-intensity focused ultrasound therapy for localized prostate cancer: An interim analysis of the multinational FASST study. *Eur J Clin Invest* 2024; 54: e14192.
56. Mala KS, Plage H, Modl L et al. Follow-Up of Men Who Have Undergone Focal Therapy for Prostate Cancer with HIFU-A Real-World Experience. *J Clin Med* 2023; 12.
57. Mortezaei A, Krauter J, Gu A et al. Extensive Histological Sampling following Focal Therapy of Clinically Significant Prostate Cancer with High Intensity Focused Ultrasound. *J Urol* 2019; 202: 717-24.
58. Nahar B, Ajami T, Williams A et al. Survival Outcomes and Recurrence Patterns Following Focal High-intensity Focused Ultrasound Treatment for Localized Prostate Cancer: Insights on Patient Selection and Lessons Learned. *Eur Urol Focus* 2025; 11: 235-41.
59. Nahar B, Bhat A, Reis IM et al. Prospective Evaluation of Focal High Intensity Focused Ultrasound for Localized Prostate Cancer. *J Urol* 2020; 204: 483-9.
60. Nyk L, Michalak W, Szeplinski S et al. High-Intensity Focused-Ultrasound Focal Therapy Versus Laparoscopic Radical Prostatectomy: A Comparison of Oncological and Functional Outcomes in Low- and Intermediate-Risk Prostate Cancer Patients. *J Pers Med* 2022; 12.
61. Reddy D, Peters M, Shah TT et al. Cancer Control Outcomes Following Focal Therapy Using High-intensity Focused Ultrasound in 1379 Men with Nonmetastatic Prostate Cancer: A Multi-institute 15-year Experience. *Eur Urol* 2022; 81: 407-13.
62. Rischmann P, Gelet A, Riche B et al. Focal High Intensity Focused Ultrasound of Unilateral Localized Prostate Cancer: A Prospective Multicentric Hemiblation Study of 111 Patients. *Eur Urol* 2017; 71: 267-73.
63. Rompre-Brodeur A, Marcq G, Tholomier C et al. Role of Systematic Control Biopsies following Partial Gland Ablation with High-Intensity Focused Ultrasound for Clinically Significant Prostate Cancer. *J Urol* 2021; 206: 1177-83.
64. Rosenhammer B, Niessen C, Rotzinger L et al. Oncological Outcome and Value of Postoperative Magnetic Resonance Imaging after Focal High-Intensity Focused Ultrasound Therapy for Prostate Cancer. *Urol Int* 2019; 103: 270-8.
65. Rosta G, Turba S, Mun DH et al. Focal Therapy Using High-Intensity Focused Ultrasound for Low- and Intermediate-Risk Prostate Cancer: Results from a Prospective, Multicenter Feasibility Trial. *Cancers (Basel)* 2025; 17.
66. S DEL, Checcucci E, Piramide F et al. MRI/real-time ultrasound image fusion guided high-intensity focused ultrasound: a prospective comparative and functional analysis of different ablative techniques. *Minerva Urol Nephrol* 2023; 75: 172-9.
67. Schmid FA, Schindele D, Mortezaei A et al. Prospective multicentre study using high intensity focused ultrasound (HIFU) for the focal treatment of prostate cancer: Safety outcomes and complications. *Urol Oncol* 2020; 38: 225-30.
68. Shoji S, Hiraiwa S, Uemura K et al. Focal therapy with high-intensity focused ultrasound for the localized prostate cancer for Asian based on the localization with MRI-TRUS fusion image-guided transperineal biopsy and 12-core transperineal systematic biopsy: prospective analysis of oncological and functional outcomes. *Int J Clin Oncol* 2020; 25: 1844-53.
69. Shoji S, Naruse J, Ohno S et al. Focal therapy using high-intensity focused ultrasound with intraoperative prostate compression for patients with localized prostate cancer: a multi-center prospective study with 7 year experience. *Prostate Cancer Prostatic Dis* 2025; 28: 782-8.
70. Sivaraman A, Marra G, Stabile A et al. Does mpMRI guidance improve HIFU partial gland ablation compared to conventional ultrasound guidance? Early functional outcomes and complications from a single center. *Int Braz J Urol* 2020; 46: 984-92.
71. Tourinho-Barbosa RR, Sanchez-Salas R, Claros OR et al. Focal Therapy for Localized Prostate Cancer with Either High Intensity Focused Ultrasound or Cryoablation: A Single Institution Experience. *J Urol* 2020; 203: 320-30.
72. van Velthoven R, Aoun F, Marcellis Q et al. A prospective clinical trial of HIFU hemiablation for clinically localized prostate cancer. *Prostate Cancer Prostatic Dis* 2016; 19: 79-83.
73. von Hardenberg J, Westhoff N, Baumunk D et al. Prostate cancer treatment by the latest focal HIFU device with MRI/TRUS-fusion control biopsies: A prospective evaluation. *Urol Oncol* 2018; 36: 401e1-e9.
74. Westhoff N, Ernst R, Kowalewski KF et al. Medium-term Oncological Efficacy and Patient-reported Outcomes After Focal High-intensity Focused Ultrasound: The FOXPRO Trial. *Eur Urol Focus* 2023; 9: 283-90.
75. Yap T, Ahmed HU, Hindley RG et al. The Effects of Focal Therapy for Prostate Cancer on Sexual Function: A Combined Analysis of Three Prospective Trials. *Eur Urol* 2016; 69: 844-51.
76. Yee CH, Chiu PK, Teoh JY et al. High-Intensity Focused Ultrasound (HIFU) Focal Therapy for Localized Prostate Cancer with MRI-US Fusion Platform. *Adv Urol* 2021; 2021: 7157973.
77. Blazeviski A, Amin A, Scheltema MJ et al. Focal ablation of apical prostate cancer lesions with irreversible electroporation (IRE). *World J Urol* 2021; 39: 1107-14.

78. Blazevski A, Scheltema MJ, Yuen B et al. Oncological and Quality-of-life Outcomes Following Focal Irreversible Electroporation as Primary Treatment for Localised Prostate Cancer: A Biopsy-monitored Prospective Cohort. *Eur Urol Oncol* 2020; 3: 283-90.
79. Collettini F, Enders J, Stephan C et al. Image-guided Irreversible Electroporation of Localized Prostate Cancer: Functional and Oncologic Outcomes. *Radiology* 2019; 292: 250-7.
80. Enikeev D, Taratkin M, Morozov A et al. Focal irreversible electroporation for localized prostate cancer management: prospective assessment of efficacy and safety. *Minerva Urol Nefrol* 2020; 72: 644-5.
81. George AK, Miodinovic R, Patel AR et al. Irreversible Electroporation for Prostate Tissue Ablation in Patients with Intermediate-risk Prostate Cancer: Results from the PRESERVE Trial. *Eur Urol* 2025.
82. Gielchinsky I, Lev-Cohain N. Focal Irreversible Electroporation for Localized Prostate Cancer - Oncological and Safety Outcomes Using mpMRI and Transperineal Biopsy Follow-Up. *Res Rep Urol* 2023; 15: 27-35.
83. Giganti F, Stabile A, Giona S et al. Prostate cancer treated with irreversible electroporation: MRI-based volumetric analysis and oncological outcome. *Magn Reson Imaging* 2019; 58: 143-7.
84. Guenther E, Klein N, Zapf S et al. Prostate cancer treatment with Irreversible Electroporation (IRE): Safety, efficacy and clinical experience in 471 treatments. *PLoS One* 2019; 14: e0215093.
85. Miñana López B, Andrés Boville G, Barbas Bernardos G et al. Focal Therapy of Prostate Cancer Index Lesion With Irreversible Electroporation. A Prospective Study With a Median Follow-up of 3 Years. *J Urol* 2023; 209: 261-70.
86. Murray KS, Ehdaie B, Musser J et al. Pilot Study to Assess Safety and Clinical Outcomes of Irreversible Electroporation for Partial Gland Ablation in Men with Prostate Cancer. *J Urol* 2016; 196: 883-90.
87. Popeneciu IV, Mohr MN, Strauss A et al. Personalized Treatment Strategy in "Low-Risk Prostate Cancer Active Surveillance Candidates" Using Irreversible Electroporation: Prospective Evaluation of Feasibility, Morbidity, Functional and Oncological Outcomes. *World J Mens Health* 2024; 42: 821-9.
88. Scheltema MJ, Geboers B, Blazevski A et al. Median 5-year outcomes of primary focal irreversible electroporation for localised prostate cancer. *BJU Int* 2023; 131 Suppl 4: 6-13.
89. Shin D, Yoon CE, Kwon HJ et al. Irreversible electroporation for prostate cancer using PSMAPET-CT. *Prostate Int* 2023; 11: 40-5.
90. Suberville M, Zhang K, Woillard JB et al. Oncological and Functional Outcomes of Hemi-Ablation Versus Focal Ablation for Localized Prostate Cancer Using Irreversible Electroporation. *Cancers (Basel)* 2025; 17.
91. Ting F, Tran M, Bohm M et al. Focal irreversible electroporation for prostate cancer: functional outcomes and short-term oncological control. *Prostate Cancer Prostatic Dis* 2016; 19: 46-52.
92. Valerio M, Dickinson L, Ali A et al. Nanoknife Electroporation Ablation Trial: A Prospective Development Study Investigating Focal Irreversible Electroporation for Localized Prostate Cancer. *J Urol* 2017; 197: 647-54.
93. van den Bos W, Scheltema MJ, Siriwardana AR et al. Focal irreversible electroporation as primary treatment for localized prostate cancer. *BJU Int* 2018; 121: 716-24.
94. Wang H, Xue W, Yan W et al. Extended Focal Ablation of Localized Prostate Cancer With High-Frequency Irreversible Electroporation: A Nonrandomized Controlled Trial. *JAMA Surg* 2022; 157: 693-700.
95. Yaxley WJ, Gianduzzo T, Kua B, Oxford R, Yaxley JW. Focal therapy for prostate cancer with irreversible electroporation: Oncological and functional results of a single institution study. *Investig Clin Urol* 2022; 63: 285-93.
96. Zhang K, Teoh J, Laguna Pet al. Effect of Focal vs Extended Irreversible Electroporation for the Ablation of Localized Low- or Intermediate-Risk Prostate Cancer on Early Oncological Control: A Randomized Clinical Trial. *JAMA Surg* 2023; 158: 343-9.
97. Zhang K, Teoh J, Laguna Pet al. The diagnostic value of MRI for persistent prostate cancer following irreversible electroporation focal therapy. *BJU Int* 2025; 136: 120-7.
98. Aker MN, Brisbane WG, Kwan L et al. Cryotherapy for partial gland ablation of prostate cancer: Oncologic and safety outcomes. *Cancer Med* 2023; 12: 9351-62.
99. Bahar P, Oh E, Desai K et al. MRI Can Replace Early Protocol Biopsy in Patients Undergoing Focal Cryoablation for Prostate Cancer. *J Endourol* 2025; 39: 725-30.
100. Baskin A, Charondo LB, Balakrishnan A et al. Medium Term Outcomes of Focal Cryoablation for Intermediate and High Risk Prostate Cancer: MRI and PSA are Not Predictive of Residual or Recurrent Disease. *Urol Oncol* 2022; 40: 451e15- e20.
101. Basourakos SP, Al Hussein Al Awamli B, Bianco FJ et al. Feasibility of in-office MRI-targeted partial gland cryoablation for prostate cancer: an IDEAL stage 2A study. *BMJ Surg Interv Health Technol* 2020; 2: e000056.
102. Bianco FJ, Martinez-Salamanca JI, Lopez Prieto A et al. Natural History of Patients Undergoing Transperineal Fusion Biopsies to Transperineal Fusion Cryoablation: 5-Year Outcomes. *J Endourol* 2025; 39: S3-s12.
103. Bossier R, Sanguedolce F, Territo A et al. Whole and hemi-gland cryoablation for primary localized prostate cancer: Short and medium-term oncological and functional outcomes. *Actas Urol Esp (Engl Ed)* 2020; 44: 172-8.
104. Fernandez-Pascual E, Manfredi C, Martin C et al. mpMRI-US Fusion-Guided Targeted Cryotherapy in Patients with Primary Localized Prostate Cancer: A Prospective Analysis of Oncological and Functional Outcomes. *Cancers (Basel)* 2022; 14.
105. Khan A, Khan AU, Siref L, Feloney M. Focal Cryoablation of the Prostate: Primary Treatment in 163 Patients With Localized Prostate Cancer. *Cureus* 2023; 15: e37172.
106. Lepor H, Cao T, Tafa M, Wysock JS. Functional, oncological, regret and complications following partial gland cryo-ablation for low-risk prostate cancer associated with mpMRI targets. *Urol Oncol* 2025; 43: 398.e1- e6.
107. Lepor H, Parry S, Tafa M, Rapoport E, Wysock JS. MRI performance to identify clinically significant prostate cancer following focal cryo-ablation. *BJUI Compass* 2025; 6: e70120.
108. Lepor H, Thakker S, Rapoport E et al. Factors predicting sexual function dissatisfaction following primary partial gland cryoablation. *J Sex Med* 2026; 23.
109. Marquis A, Schwartz MJ, Armache A et al. Novel MRI-US fusion with advanced annotation in focal cryoablation for prostate cancer: medium-term outcomes and learning curve insights. *BJU Int* 2026; 137: 189-200.
110. Selvaggio O, Finati M, Falagario UG et al. Treatment of localized prostate cancer in elderly patients: the role of partial cryoablation. *Int Urol Nephrol* 2023; 55: 1125-32.
111. Shah TT, Peters M, Eldred-Evans D et al. Early-Medium-Term Outcomes of Primary Focal Cryotherapy to Treat Nonmetastatic Clinically Significant Prostate Cancer from a Prospective Multicentre Registry. *Eur Urol* 2019; 76: 98-105.
112. Sidana A, Tayebi S, Blank F et al. Magnetic resonance imaging-ultrasound fusion guided focal cryoablation for men with intermediate-risk prostate cancer. *Urol Oncol* 2024; 42: 158e1- e10.
113. Tan WP, Chang A, Sze C, Polascik TJ. Oncological and Functional Outcomes of Patients Undergoing Individualized Partial Gland Cryoablation of the Prostate: A Single-Institution Experience. *J Endourol* 2021; 35: 1290-9.
114. Tan YG, Law YM, Ngo NT et al. Patient-reported functional outcomes and oncological control after primary focal cryotherapy for clinically significant prostate cancer: A Phase II mandatory biopsy-monitored study. *Prostate* 2023; 83: 781-91.
115. Tay KJ, Hong BH, Ong EHW et al. Transcriptomic predictors of prostate cancer recurrence following focal cryotherapy: a pooled analysis of phase II trial and prospective cohort data. *J Natl Cancer Cent* 2025; 5: 515-23.
116. Wysock JS, Rapoport E, Hernandez H, Gogaj R, Lepor H. Biopsy Assessment of Oncologic Control 3 Years Following Primary Partial Gland Cryoablation: A Prospective Cohort Study of Men With Intermediate-risk Prostate Cancer. *J Urol* 2023; 210: 454-64.
117. Light A, Peters M, Arya M et al. Salvage Focal Therapy vs Radical Prostatectomy for Localized Radiorecurrent Prostate Cancer. *JAMA Oncol* 2026; 12: 364-73.

References: Table 2

1. Wei S, Li C, Li M, et al. Radioactive Iodine-125 in Tumor Therapy: Advances and Future Directions. *Front Oncol*. 2021 Sep 30;11:717180.
2. Shehata Elhelfi A, Albahar H, Shah U, et al. High intensity focused ultrasound: The fundamentals, clinical applications and research trends. *Diagnostic and Interventional Imaging*. 2018;99: 349-359.
3. Geboers B, Scheltema MJ, Jung J et al. Irreversible electroporation of localised prostate cancer downregulates immune suppression and induces systemic anti-tumour T-cell activation – IRE-IMMUNO study. *BJU Int*. 2025; 135:319-328
4. Baust J G, Gage A A, Bjerklund Johansen T E, et al. Mechanisms of cryoablation: clinical consequences on malignant tumors. *Cryobiology*. 2014;68: 1-11.
5. Adhami M, Cheng J, Anderson E et al. Functional outcomes and toxicity following focal low-dose-rate brachytherapy for low-intermediate risk prostate cancer. *BJUJ Compass* 2026; 7: e70119.
6. Durazo-Ruiz FO, Roldán-Testillano R, Rodríguez-Sanchez L, et al. Irreversible electroporation in focal therapy for prostate cancer: current status and future directions. *Prostate International*. 7 November 2025
7. Henry A, Pieters B R, Siebert F-A et al. GEC-ESTROACROP prostate brachytherapy guidelines. *Radiotherapy and Oncology* 1672022; 244–251.
8. Ranjan P, Saurabh G, Bansal R, et al. High-Intensity Focused Ultrasound vs. Cryotherapy as Primary Treatment for Prostate Cancer. *Indian J Urol*. 2008; 24: 16–21.
9. Narayanan G. Irreversible Electroporation. *Semin Intervent Radiol*. 2015;32: 349–355.
10. Lee EW, Huang W C. Minimally Invasive Ablative Therapies for Definitive Treatment of Localized Prostate Cancer in the Primary Setting. *Prostate Cancer*. 2011; :394182.
11. Stone N N, Stock R G. Prostate brachytherapy in men with gland volume of 100 cc or greater: Technique, cancer control, and morbidity. *Brachytherapy* 2013; 12: 217-221.
12. Yu Z, Wen Q, Guo Y, et al. Preliminary clinical outcomes of irreversible electroporation in the management of prostatic hyperplasia. *Medicine (Baltimore)*. 2025; 28; :e46166
13. Roberts C B, Jang T L, Shao Y-H, et al. Treatment profile and complications associated with cryotherapy for localized prostate cancer: A population-based study. *Prostate Cancer Prostatic Dis*. 2011;14:313–319.
14. van den Bos W, de Bruin D M, Muller B G, et al. The safety and efficacy of irreversible electroporation for the ablation of prostate cancer: a multicentre prospective human in vivo pilot study protocol. *BMJ Open*. 2014;4:e006382.
15. NICE INTERVENTIONAL PROCEDURES PROGRAMME. Interventional procedures overview of high-intensity focused ultrasound for prostate cancer. <https://www.nice.org.uk/guidance/htg70/evidence/overview-pdf-304022125>
16. Faiella E, Santucci D, Vertulli D, et al. Irreversible Electroporation (IRE) for Prostate Cancer (PCa) Treatment: The State of the Art. *J Pers Med*. 2024; 14:137.
17. Paesano N, Quintas J J, Barreras S G, et al. High-intensity focused ultrasound (HIFU) for prostate cancer - step by step technique. *Urology Video Journal*. 2026; 29: 100383.
18. Cohen J K, Miller R J, Shuman B A. Urethral warming catheter for use during cryoablation of the prostate. *Urology*, 1995;45: 861-864.
19. Cytron S, Paz A, Kravchick S, et al. Active Rectal Wall Protection Using Direct Transperineal Cryo-Needles for Histologically Proven Prostate Adenocarcinomas. *European Urology*, 2003; 44: 315-321.
20. Levy D A. Correlation of Thermocouple Data With Voiding Function After Prostate Cryoablation. *Urology*, 2010; 75: 482-486.
21. Yee C-H, Chiu P K-F, Teoh J Y-C, et al. High-Intensity Focused Ultrasound (HIFU) Focal Therapy for Localized Prostate Cancer with MRI-US Fusion Platform. *Adv Urol*. 2021; :7157973.
22. Baur A D J, Colletini F, Enders J, et al. MRI-TRUS fusion for electrode positioning during irreversible electroporation for treatment of prostate cancer. *Diagn Interv Radiol*. 201;23: 321–325.
23. Fernández-Pascual E, Manfredi C, Martín C, et al. mpMRI-US Fusion-Guided Targeted Cryotherapy in Patients with Primary Localized Prostate Cancer: A Prospective Analysis of Oncological and Functional Outcomes. *Cancers (Basel)* 2022; 14:2988.
24. Anderson E, Smyth L M L, O’Sullivan R, et al. Focal low dose-rate brachytherapy for low to intermediate risk prostate cancer: preliminary experience at an Australian institution. *Transl Androl Urol*. 2021;10: 3591–3603.
25. Wendler J J, Fischbach K, Rieke J, et al. Irreversible Electroporation (IRE): Standardization of Terminology and Reporting Criteria for Analysis and Comparison. *Pol J Radiol*, 2016; 81:54-64.
26. Lodeizen O, de Bruin M, Eggen S, et al. Ablation energies for focal treatment of prostate cancer. *World J Urol*. 2018; 37: 409–418.
27. Tan Y G, Law Y M, Ngo N T, et al. Patient reported functional outcomes and oncological control after primary focal cryotherapy for clinically significant prostate cancer: A Phase II mandatory biopsy-monitored study. *Prostate*. 2023; 83: 781–791.
28. Brisbane W G, Priester A M, Ballon J, et al. Targeted prostate biopsy: umbra, penumbra, and value of perilesional sampling. *Eur Urol*. 2022; 82: 303–310.
29. Lantz A, Nordlund P, Falagarío U, et al. Prostate Cancer IRE Study (PRIS): A Randomized Controlled Trial Comparing Focal Therapy to Radical Treatment in Localized Prostate Cancer. *Eur Urol Open Sci*. 2023; 51: 89–94.
30. Nassiri N, Chang E, Lieu P, et al. Focal Therapy Eligibility Determined by Magnetic Resonance Imaging/ Ultrasound Fusion Biopsy. *J Urol* 2018;199:453-8. 10.1016/j.juro.2017.08.085
31. Sze C, Tsivian E, Tay K J, et al. Anterior gland focal cryoablation: proof-of-concept primary prostate cancer treatment in select men with localized anterior cancers detected by multi-parametric magnetic resonance imaging. *BMC Urol*. 2019; 5: 19:127.
32. Rakauskas A, Marra G, Heidegger I, et al. Focal Therapy for Prostate Cancer: Complications and Their Treatment. *Front Surg*. 2021; 8: 696242.
33. de Castro Abreu A L, Ashrafi A N, Gill I S, et al. Contrast-Enhanced Transrectal Ultrasound for Follow-up After Focal HIFU Ablation for Prostate Cancer. *J Ultrasound Med*. 2019;38: 811-819.
34. Colletini F, Enders J, Stephan C, et al. Image-guided Irreversible Electroporation of Localized Prostate Cancer: Functional and Oncologic Outcomes. *Radiology*. 2019;292: 250-257.



As a trusted partner in prostate cancer care, BXTA delivers innovation in cancer diagnosis pathway and provides proven low dose rate (LDR) brachytherapy products and services. Our mission is to empower healthcare professionals worldwide to treat prostate cancer, help extend lives and improve the quality of outcomes for patients.

Please contact your local BXTA celyon representative for the most current information or contact us on the address listed below.

International Headquarters

BXTA
The Coach House, Grenville Court
Britwell Road, Burnham, SL18DF, United Kingdom
Phone: +44 01628 308 420

info@bxta.com

All images and photographs © 2026 BXTA celyon Limited and its subsidiaries.

BXTA, the BXTA Logo are trademarks of BXTA celyon Ltd. All other company brand, product, and service names may be trademarks or registered trademarks of their respective holders.

All rights reserved.