



**For smaller patient access,
keep it short and sweet.**

Created by and for neonatal and pediatric clinicians, the MicroSlide® family offers appropriately sized products/tools for smaller patients. The MicroSlide® pediatric MST product line offers 2F and 3F sheaths in a variety of lengths and configurations. Offered as a full kit, or as a separate sterile component, the sheaths can be paired with appropriate .008", .010", .014" or .018" guidewire and 24G 2.5cm or 4cm echogenic needle.



Theragenics®



theragenics.com



info@theragenicsinterventional.com



800.639.2800

Galt Medical Corp. is a subsidiary of Theragenics Corporation

Theragenics Corporation, 5203 Bristol Industrial Way, Buford, Georgia 30518, USA

MKT-100-35 Rev. E



MicroSlide has been a life changer for our pediatric population! It has helped eliminate multiple sticks greatly.

*Katherine Sendelback,
PICC, LeBonheur & Mid South*



Excellent Product

IR, MD



Tearaway Introducer Kits

Catalog No.	Sheath Size	Guidewire Size	Guidewire Material	Echo Needle	Units/Box
KIT-051-00	2F 2F x 3cm	.008" x 20cm	NT/NT	24G x 2.5cm	5
KIT-051-01	2F 2F x 3cm	.008" x 40cm	NT/NT	24G x 2.5cm	5
KIT-051-02	2F 2F x 3cm	.008" x 20cm	NT/NT	24G x 4cm	5
KIT-051-03	2F 2F x 3cm	.008" x 40cm	NT/NT	24G x 4cm	5
KIT-051-04	3F 3F x 3cm	.010" x 20cm	NT/NT	24G x 2.5cm	5
KIT-051-05	3F 3F x 3cm	.010" x 40cm	NT/NT	24G x 2.5cm	5
KIT-051-06	3F 3F x 3cm	.010" x 20cm	NT/NT	24G x 4cm	5
KIT-051-11	3F 3F x 3cm	.010" x 40cm	NT/NT	24G x 4cm	5
KIT-051-13	2F 2F x 3cm	.014" x 20cm	NT/Au	24G x 4cm	5
KIT-051-14	2F 2F x 3cm	.014" x 40cm	NT/Au	24G x 4cm	5

Tearaway Introducer Sets

Catalog No.	Sheath Size	Accepts Guidewire	Units/Box
INT-101-02	2F 2F x 3cm	.008"	10
INT-101-14	2F 2F x 3cm	.014"	10
INT-101-15	2F 2F x 7cm	.014"	10
INT-101-03	3F 3F x 3cm	.018"	10

Set Contents: PTFE Sheath and Dilator only

Guidewires

Catalog No.	Size	Guidewire Material	Units/Box
SGW-051-06	.008" x 20cm	NT/NT	5
SGW-051-08	.008" x 40cm	NT/NT	5
SGW-051-03	.008" x 60cm	NT/NT	10
SGW-051-07	.010" x 20cm	NT/NT	5
SGW-051-09	.010" x 40cm	NT/NT	5
SGW-051-04	.010" x 60cm	NT/NT	10
SGW-056-62	.014" x 20cm	NT/Au	5
SGW-056-60	.014" x 40cm	NT/Au	5

Needles

Catalog No.	Size	Length	Tip	Units/Box
NDL-100-29	24G	2.5cm	echogenic	25
NDL-100-30	24G	4cm	echogenic	25

