

Brachytherapy solutions to better serve physicians and their patients

## TheraSeed<sup>®</sup> & AgX100<sup>®</sup>

### Pre-loaded Needle Kit options to meet your LDR brachytherapy implant needs...

#### Customized per a physician's plan

- Manufactured to a patient's individual treatment plan
- Stranded seed configurations
  - TheraStrand<sup>®</sup>
  - TheraSleeve<sup>®</sup>
- Non-stranded seed configurations
  - TheraLoad<sup>®</sup>
  - AnchorLoad<sup>®</sup>
- Variable seed and spacer configurations
  - 2.75 mm, 5 mm & 1.0 cm spacer lengths
  - Stranded: 6.5 cm & 7.0 cm maximum configured length
  - Non-stranded: 7.0 cm maximum configured length
- Round or square hub 18-gauge implant needle options
  - Marked with grid locations
  - Wax or synthetic-tension plug
  - Stylet lock prevents accidental seed deployment
- Vertical or horizontal needle dispensing positions



#### Quality assured

- Sterile & shielded
- Patient specific information label
- Needle configuration plan
- Verification radiograph documents exact seed position per needle
- Manufacturers Certificate of Calibration; an additional independent seed assay available upon request



## Pre-loaded Needle Kit packaging



## Needle Trays

- Sterile
- Radio-protective shielding
- Qty 1-16 needles per tray
- Qty 1-3 needle trays per shipper



## Shipping Container

- Distinctive Theragenics® blue & gold colored box
  - Hassle-free shipment identification
  - Segregated section accommodates loose seed and seed cartridge shielded containers as needed
  - Dimensions: 61cm x 32cm x 20cm (24in x 13in x 8in)



**For specific product information and Instructions for Use, please consult your Theragenics Brachytherapy Specialist or Customer Service.**

### Customer Service

Phone: 877-444-7333

Fax: 800-458-4303

Email: [customerservice@theragenics.com](mailto:customerservice@theragenics.com)

Hours of Operation: Mon - Fri, 8am - 6pm ET

### Theragenics Corporation

5203 Bristol Industrial Way

Buford, Georgia 30518

USA