

Precision and Performance

The ELITE HV® R transitionless radial introducer sheath is designed for precise access, smooth entry, and consistent performance. Featuring transitionless, hydrophilic sheath properties, combined with a premium guidewire and needle, the ELITE HV® R kit provides the tools needed for successful, atraumatic radial access.

 **23% Less Insertion Force**
than a comparable sheath*

 **Hydrophilic Coating**
for easy insertion and removal


Transitionless Sheath Tip
for smooth entry


*1:1 simulated data on file

Radial Artery Access Sheath

Sheath/dilator assembly accommodates an .018" guidewire and is available in sizes from 4F to 7F and in 7cm and 11cm lengths.




Sheath Tip
ID/Sheath OD
Specifications



Guidewire

Kits include an .018" x 45cm guidewire with nitinol core and gold radiopaque tip (shown below) or stainless tip and core. Separate sterile guidewires in additional lengths are also available.



Catalog No.	Diameter	Length	Description	Box Qty.
SGW-056-01	.018"	40cm	NT/Au	5
SGW-056-04	.018"	45cm	NT/Au	5
SGW-015-02	.018"	40cm	SS/SS	10
SGW-015-08	.018"	45cm	SS/SS	10

21G Needle

Kits include a premium 21G x 2.5cm needle, designed for precision and featuring a bevel indicator on the hub. Separate sterile needles, including X Sharp[®], are also available.



Catalog No.	Needle Size	Description	Box Qty.
NDL-100-32	21G 21G x 2.5cm	Smooth	25
NDL-100-17	21G 21G x 4cm	Smooth	25
NDL-107-14	21G 21G x 4cm	Smooth, X Sharp	25

Catalog No.	Sheath Size	Dilator Length	Wire	Needle Size	Sheath ID	Sheath OD	Tip ID
CLI-230-45	4F 4F x 7cm	11.8"	.018" x 45cm NT/Au	21G 21G x 2.5cm	0.066"	0.089"	0.061"
CLI-230-55	5F 5F x 7cm	11.8"	.018" x 45cm NT/Au	21G 21G x 2.5cm	0.079"	0.102"	0.074"
CLI-230-06	6F 6F x 7cm	11.8"	.018" x 45cm NT/Au	21G 21G x 2.5cm	0.086"	0.109"	0.086"
CLI-230-07	7F 7F x 7cm	12.0"	.018" x 45cm NT/Au	21G 21G x 2.5cm	0.099"	0.122"	0.094"
CLI-228-45	4F 4F x 7cm	11.8"	.018" x 45cm SS/SS	21G 21G x 2.5cm	0.066"	0.089"	0.061"
CLI-228-55	5F 5F x 7cm	11.8"	.018" x 45cm SS/SS	21G 21G x 2.5cm	0.079"	0.102"	0.074"
CLI-228-06	6F 6F x 7cm	11.8"	.018" x 45cm SS/SS	21G 21G x 2.5cm	0.086"	0.109"	0.086"
CLI-228-07	7F 7F x 7cm	12.0"	.018" x 45cm SS/SS	21G 21G x 2.5cm	0.099"	0.122"	0.094"
CLI-231-45	4F 4F x 11cm	15.8"	.018" x 45cm NT/Au	21G 21G x 2.5cm	0.066"	0.089"	0.061"
CLI-231-55	5F 5F x 11cm	15.8"	.018" x 45cm NT/Au	21G 21G x 2.5cm	0.079"	0.102"	0.074"
CLI-231-06	6F 6F x 11cm	15.8"	.018" x 45cm NT/Au	21G 21G x 2.5cm	0.086"	0.109"	0.086"
CLI-231-07	7F 7F x 11cm	16.0"	.018" x 45cm NT/Au	21G 21G x 2.5cm	0.099"	0.122"	0.094"
CLI-229-45	4F 4F x 11cm	15.8"	.018" x 45cm SS/SS	21G 21G x 2.5cm	0.066"	0.089"	0.061"
CLI-229-55	5F 5F x 11cm	15.8"	.018" x 45cm SS/SS	21G 21G x 2.5cm	0.079"	0.102"	0.074"
CLI-229-06	6F 6F x 11cm	15.8"	.018" x 45cm SS/SS	21G 21G x 2.5cm	0.086"	0.109"	0.086"
CLI-229-07	7F 7F x 11cm	16.0"	.018" x 45cm SS/SS	21G 21G x 2.5cm	0.099"	0.122"	0.094"

Quantity: 5 per box

