

Fiducial Markers for Radiation Therapy

AnchorMarker[®]
TECHNOLOGY



Markers are not to scale or size.

Radiopaque Gold markers for identification of anatomic regions by providing reference position for a proposed treatment site.

AnchorMarker[®] Needles

Gold markers encapsulated with Theragenics' proprietary, synthetic bio-absorbable polymer are designed to:

- Improve fixity in tissue¹
- Reduce misalignment¹
- Minimize migration and displacement²
 - 0.8mm diameter x 4.5mm length gold marker
 - Sterile, pre-loaded 18G x 20cm needle
 - Available in packs of 1 or 3 needles

Gold Fiducial Marker Needles

Packaging Options:

1.0mm diameter Gold Marker

- 3mm or 5mm lengths
- Sterile, pre-loaded 18G x 20cm needle
- Available in packs of 1 or 3 needles

1.2mm diameter Gold Marker

- 3mm or 5mm lengths
- Sterile, pre-loaded 17G x 20cm needle
- Available in packs of 1 or 3 needles

References

1. Badwan H, Shanahan A, Adams M, et al. AnchorSeed for the reduction of source movement in prostate brachytherapy with the Mick™ Applicator implant technique, *Brachytherapy* 9 (2010) 23-26
2. Sarkar A, Donovanik V, Zhang I, et al. Prostate implant dosimetric outcomes and migration patterns between bio-absorbable coated and uncoated brachytherapy seeds, *Brachytherapy* 12 (2013) 356-361

Content intended for the USA. Indications may vary in different countries; please either consult the appropriate IFU or your local distributor. Theragenics Corporation and AnchorMarker are registered trademarks of Theragenics Corporation. ©2021 Theragenics Corporation

For ordering information, please consult your Theragenics Brachytherapy Specialist or Customer Service.

Theragenics



theragenics.com



customerservice@theragenics.com



877.444.7333

Theragenics Corporation, 5203 Bristol Industrial Way, Buford, Georgia 30518, USA

ref033121